Problem	Potential impact on patient / device	Presumed cause	Recommended action
An allergic reaction of the patient	-	<i>The patient is sensitive to allergenic compounds contained in the product</i>	<i>Use other material suitable for the patient</i>
	-	<i>Increased level of residual monomer</i>	<i>Strictly follow the powder/liquid ratio</i> <i>Use recommended way of reduction residual monomer level in the IFU.</i>
Porous structure in the finished prosthesis	<i>Increased residual monomer level - may cause irritation and allergic reaction</i>	<i>Incorrect powder/liquid ratio</i>	<i>Strictly follow the powder/liquid ratio</i>
		<i>Incorrect process of pressing the acrylic dough in the flask under the hydraulic press</i>	<i>Use the pressure and pressing time specified in the IFU.</i>
	<i>Reduction of the mechanical properties of the prosthesis, breaking of the prosthesis during use, teeth falling out of the denture</i>	<i>Incorrect polymerization process</i>	<i>Use the recommended process parameters.</i> <i>In the case of thick lower dentures, use a prolonged polymerization</i>
		<i>Not enough of acrylic material in the flask</i>	<i>10-15% more acrylic mass should be put in the flask. The excess material coming out of the can must be visible.</i>
Reduction of the mechanical properties of the prosthesis after polymerization	<i>Breaking dentures during use</i>	<i>Incorrect powder/liquid ratio</i>	<i>Strictly follow the powder/liquid ratio</i>
		<i>Incorrect polymerization process</i>	<i>Use the recommended process parameters.</i>
Rough surface of the prosthesis on the mucosal surface	<i>Increased affinity for the accumulation of fungi and microbial flora - possible irritation and microbial infections</i>	<i>Incorrect isolation of stone model</i>	<i>Check the type of isolator used. A film-forming insulator, eg Izo-sol, should be used.</i>
		<i>Inaccurate scalding process of the wax denture</i>	<i>During the scalding of the wax denture, do not allow the wax to liquefy completely, which will be in result absorbed into the gypsum structure and preventing its proper isolation with a film-forming preparation.</i>
Incorrect denture color	<i>The esthetics of the prosthesis is not acceptable</i>	<i>Incorrect powder/liquid ratio</i>	<i>Strictly follow the powder/liquid ratio</i>
Change in color and discoloration of the denture during use	<i>The esthetics of the prosthesis is not acceptable</i>	<i>Soaking the technician's hands with water / using moisturizing creams when processing the material</i>	<i>Use polyethylene HDPE gloves</i>
		<i>Usage of aggressive agents for cleaning dentures containing oxidizing substances</i>	<i>Advise the patient to use other, less aggressive agents specified in the IFU.</i>
		<i>Sludge from smoking, drinking coffee, dyes contained in food or spices used in food, medications</i>	<i>Use of denture cleaning agents recommended by the IFU</i>

Falling out acrylic teeth after polymerization	<i>Damage of the prosthesis. Loss of intended properties and functions.</i>	<i>Leaving wax residues on the surface of the teeth</i>	<i>Thorough removal of wax residues from acrylic teeth after scalding</i>
		<i>Incorrect preparation of the teeth for connection with the denture.</i>	<i>Performing retention on the mucosal surface of an acrylic tooth Soaking the prepared surfaces of acrylic teeth with monomer</i>
		<i>Incorrect polymerization process</i>	<i>Use the recommended process parameters.</i>
The denture plate does not fit onto stone model	<i>Denture retention problems</i>	<i>Polymerization shrinkage</i>	<i>Strictly follow the powder/liquid ratio Use the recommended process parameters.</i>
Incorrect occlusion		<i>Incorrect consistency of the acrylic (too thick) when putting into the flask</i>	<i>Reducing the waiting time after mixing powder and liquid Lowering the temperature in the laboratory</i>
		<i>Incorrect process of pressing the acrylic dough in the flask under the hydraulic press</i>	<i>Use the pressure and pressing time specified in the IFU. Use clean and sealed polymerization flask</i>

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