## VILLACRYL STC HOT INSTRUCTIONS FOR USE - EN

**Heat-curing acrylic resin for crowns and bridges veneering**

**Product characteristics:** heat-cured acrylic resin for crowns and bridges facing.

**Classification accordance EN ISO 10477:** Type 2, Class 1

**Indications for use:** acrylic faced crowns and bridges.

**Clinical benefits**

Restoration of biting and chewing functions lost as a result of partial or complete loss of teeth.

**Target group or target groups of patients and intended users of the devices**

**Users:** Professionals in the dental technology industry - certified dental technicians producing crowns and bridges.

**Patients:** People of different ages and sexes, with partial or complete lack of teeth, eligible for use of acrylic crowns or bridges

**Contraindications:** Do not use the product at persons sensitized to any of the ingredients.

**Do not use** the product by pregnant women and nursing mothers in the group of dental technicians.

**Side effects:** Possible irritation and inflammation of the oral cavity caused by fungi and bacteria adhering to the crown or bridge in the absence of appropriate hygiene.

**Dental technicians**

Irritation, redness or symptoms of hypersensitivity, including the occurrence of anaphylactic shock, in case of allergy to any of the product ingredients.

**Patients**

Irritation, inflammation and allergic reactions caused by residual monomer washed out of the crown.

**Data necessary to identify the product and the contents of the package**

Reference number, product name, color, package size, batch number and expiry date - imprint on the product packaging.

Package content - see table below.

**Colour shades**

**Dentine:** A1; A2; A3; A3,5; A4; B1; B2; C2; C4; D2

**Enamel:** E1; E2; E3; E4

Asortment	
REF	Set
V210A1Z01	VILLACRYL STC HOT kolor A1 (80 g powder + 40 ml liquid)
V210A2Z02	VILLACRYL STC HOT kolor A2 (80 g powder + 40 ml liquid)
V210A3Z03	VILLACRYL STC HOT kolor A3 (80 g powder + 40 ml liquid)
V210A35Z04	VILLACRYL STC HOT kolor A3,5 (80 g powder + 40 ml liquid)
V210A4Z05	VILLACRYL STC HOT kolor A4 (80 g powder + 40 ml liquid)
V210B1Z06	VILLACRYL STC HOT kolor B1 (80 g powder + 40 ml liquid)
V210B2Z07	VILLACRYL STC HOT kolor B2 (80 g powder + 40 ml liquid)
V210C2Z08	VILLACRYL STC HOT kolor C2 (80 g powder + 40 ml liquid)
V210C4Z09	VILLACRYL STC HOT kolor C4 (80 g powder + 40 ml liquid)
V210D2Z12	VILLACRYL STC HOT kolor D2 (80 g powder + 40 ml liquid)
V210Z13	VILLACRYL STC HOT dentine powder 40g color A2, A3, A3,5, A4, B2, C2 + dentine powder 20g color A1, B1, C4, D2 + enamel powder 20g color E1, E2, E3, E4 + liquid 200ml
	<b>Liquid</b>
V210L02	VILLACRYL STC HOT Liquid 200ml

**Chemical composition**

**Powder:** acrylic polymers, benzoyl peroxide, pigments

**Liquid:** methyl methacrylate, cross-linking agent

**Compatible products**

For use with metals and dental alloys.

**Method of use**

**Attention:** In the purpose of the proper selection of tooth shade color it is recommended using the table attached to the instruction. The table is only a proposal of application color blends. Color should be selected individually and in dubious cases the own sample should be make.

**Preliminary works**

The metallic structure of crown or bridge shall be made by traditional method. The chewing surfaces shall be left metallic. The part designed for facing shall have retention elements (e.g. balls). When structure is polished, designs designed for facing shall be sand blasted with the jet of sand (aluminum oxide), grain size 250 μm, under pressure 4-6 bar and masked using product Villacryl Opaker constituting simultaneously the first layer of color reconstructed tooth. It is possible to mix powders of opaker creating individually first layer of color. The prepared structure place on working model and model with wax the shape of the future facing. Wax shall be replaced for acrylate by traditional method. The metal part of structure shall be embedded in dass II gypsum (Stodent II is recommended) in the lower part of polymerization box so that the faced part remains exposed.

The gypsum surface shall be isolated with gypsum–gypsum insulation after hardening. Then place the second part of the box, the so called “backing”, and it shall be also filled with class II gypsum. After 20 minutes, when gypsum is bounded to the proper hardness, polymerization box shall be opened and wax shall be scalded with boiling water. The form shall be isolated with gypsum–acrylate isolator (e.g. Izo-Sol). When the polymerization box is cooled to ambient temperature, ca 23 °C, the acrylic dough can be applied.

**Recommended following mixing proportions:** 2,4 g of power and 1 g of liquid **Dosing and mixing**

It is recommended to use gloves from polyethylene HDPE for protection of the hands.

Prepare two chemically resistant containers (glass, polyethylene, etc.) for mixing powder with liquid. Pour proper amount of monomer to the first container, add powder (dentine) of the proper color till saturation is obtained, and thoroughly mix. It is possible to mix powders of dentine creating individually the appropriate shade of color. After 2 minutes pour to the second container the proper amount of monomer, add powder (enamel) and thoroughly mix. Both containers shall be accurately covered. The acrylic paste binding level shall be controlled after short period of time.

**Application**

Paste can be put to the mould after ca 8-10 minutes from mixing powder with liquid. Thoroughly mix the dentine part, squeezing with hands protected by gloves from HDPE, place in mould and cover with polyethylene foil. Close the polymerization box and place under the press. Carefully press till the box edges are completely folded. Take the box from the press, remove the foil, cut off obliquely with sharp knife dentine part toward the secant edge. Apply the prepared previously enamel layer in place of removed acrylate, cover with foil and press the polymerization box till 3000 kg pressure. Open the box and check correctness of dentine transition into enamel. If everything is OK, it is recommended to remove the excess of the acrylate pressed aside, again assemble the box, this time without foil separator, and again press till 3000 kg pressure. Take the box from the press after 10 minutes and place in polymerization frame.

**Polymerization**

The polymerization box fastened in frame place in water with temperature 80-100°C and heat for 30 minutes.

**Processing**

The polymerization box after cooling to ambient temperature - shall be opened and the prepared gypsum element shall be released. Processing and polishing by standard method.

**Safety instructions**

**Warnings**

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

P210 Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

**Precautions**

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: The area of product contact with skin wash thoroughly with soap with water.

P501 Dispose of contents/container to an approved waste disposal plant.

Inflammable liquid (contains methyl methacrylate) - flash point +10 °C. Keep away from fire and radiation sources. While working with the material do not: consume meals or smoke tobacco. Avoid breathing vapour of liquid and direct contact of product components or the acrylic dough on the skin and eyes. Work should be carried out in good ventilated premises. Use standard ventilation. In case of contact of the powder particles with eyes the mechanical irritation may occur. In case of contact acrylic powder or liquid with eyes wash it with water and contact physician. The area of contact the product components or wet acrylic dough with skin wash thoroughly with soap water. In case of accidental swallowing contact the physician. Liquid dangerous to the environment. Do not pour into drainage system.

**Hygiene, sterilization, disinfection**

**Product**

The product is intended for multiple use. It is not supplied sterile and does not require sterilization or disinfection prior to use. To protect the product against contamination, apply the principles of good laboratory practices in dental laboratories. **The crowns made of materials** There is a risk of transmission of infection. The crown should be cleaned and disinfected by the dentist before being placed in the patient’s mouth for the first time.

For the disinfection of crown made of materials, it is recommended to use disinfectants for disinfecting dental prostheses, following the procedures recommended by the manufacturer.

It is also possible to use agents based on at least 65% ethyl or propyl alcohol. The crown should be sprayed with a disinfectant solution, and then gently wiped dry with a paper tissue. You can also use ready-to-use disposable wipes soaked in a disinfecting solution. Additionally, the crown should be washed with warm soapy water. Provide the appropriate information to the dentist.

**Information to be provided to patients / prosthesis users**

If any symptoms of an allergic reaction or irritation occur, discontinue using the crown immediately and consult a doctor.

Strong oxidizing compounds contained in cleaning agents and mouthwashes may change the color of the material after polymerization.

For the hygiene of crowns and bridges, popular and generally available toothpastes, threads, interdental brushes can be used, it is recommended to use a pressure device of the WATERPIK type. The material is not resistant to long-term effects of alcoholic beverages.

The material may change color or lose its original physical properties when exposed to chemicals or if it is misused.

Report any damage of the material to the dentist. Repairs may be only performed by a dental technician.

If a piece of damaged crown is swallowed, consult a doctor immediately. The material after polymerization is not visible on X-ray pictures.

**Storage and handling**

Liquid is sensitive to influence of temperature and UV radiation under the influence of which the uncontrolled reaction of polymerization may occur. Product components shall be stored in original containers, in dark and good ventilated compartment with temp. 5-25 °C, place not accessible by children. Do not use the product after expiry date specified on the sales packaging and labels of the product ingredients.

**Disposal considerations**

Dispose of in accordance with applicable regulations. Carton packaging, label instructions – the possibility of recycling, dispose of as communal not dangerous waste. Powder component and a packaging are not dangerous. Dispose of as communal not dangerous waste. In case of plastic container the possibility of recycling. Liquid and liquid packaging - hazardous waste. Pass to an authorized waste contractor.

**Problem prevention (2)**

Preventing problems (allergic reactions, physical and mechanical properties of the material after polymerization) - see website: **www.everall7.pl**, bookmark: products=dental acrylics → Villacryl STC Hot → downloads.

**Reporting of serious incidents:** Any serious incidents related to the products must be reported to the manufacturer and the competent authority in accordance with local regulations.

Instructions in written form, oral or information relayed during Producers presentation result from actual level of dental technique and technology. It should be treated as obligatory information. Information’s obtained from third persons (due to their different level of authorization) do not release the user from personal check whether the preparation is intended to the designated application. Due to the fact that the application takes place without the possibility of control from the producer’s side; therefore it limits the producers liability to the value of provided material.

We shall be very grateful for all comments related to products, please mail it to our address.

**MD** The symbol means a medical device as defined in EU Regulation 2017/745.

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