FORM PLAST

INSTRUCTIONS FOR USE - EN

Self-curing acrylic resin for patterns

Product characteristics: Self-curing acrylic resin for patterns. Burns without residues.

Designation

Supporting works in dental laboratories

- modeling of secondary parts of telescopic crowns, conical crowns and others,
- modeling of inserts (inlay, onlay), tongue arch, palatine arch, clasps and other elements of skeletal prosthesis,
- modeling of Maryland type adhesive bridges,
- modeling of crown-root inserts,
- structure modeling in post-implantation prosthesis,
- temporary connections in soldering.

Contraindications: In people sensitive to the ingredients of the product allergic reaction may occur. If you notice any symptoms of allergic should stop using the product and consult a doctor.

REF	Set
V220Z01	Form Plast powder 30g + liquid 2×12ml
V220Z02	Form Plast powder 100g + liquid 2×50ml
Liquid	
V220L02	Form Plast liquid 12ml

TP037	Form Plast Cup 12 pieces
TP038	Form Plast Brush

Chemical composition

Powder: acrylic polymers, benzoyl peroxide, pigments

Liquid: methyl methacrylate, ethyl methacrylate, ethylene glycol dimethacrylate, N,N-dimethyl-p-toluidne

Instructions for use

Procedure in technical and dental laboratories

Preparation: Class III and IV gypsum shall be used for working models, **Stodent III** and **Stodent IV** are particularly recommended. The gypsum model shall be isolated with compound which does not leave film on surface (DO NOT USE ISOLATORS BASED ON ALGINATES). Prepare two containers (glass, polyethylene (PE), ceramic containers or silicone cups). Transfer small amounts of powder and liquid to the prepared containers.

Procedure: Material shall be applied on model with pencil brush. Soak the tip of pencil brush in liquid and then in small amount of powder. Place acrylic drop prepared in such manner on model. The quick gelation process causes that material does not flow from the place where it has been located. Spread material on the modeled element with pencil brush. This procedure shall be repeated till the required shape and thickness of the formed element is achieved. Material is hardened automatically in ambient temperature ca 23 °C after ca 3-5 min. In case of modeling on metal (e.g. telescopic crowns), material shall be applied directly on the polished part, without any isolators. The pencil brush shall be washed in liquid, when necessary. The pencil brush shall be very carefully wash with liquid (methyl matacrylic) after completion of work.

ATTENTION: Do not soak pencil brush with powder in bottle with liquid.

The introduction of powder grains into liquid may cause not controlled polymerization process and the loss of properties by liquid.

Processing and finishing: Carefully remove the formed mass from model (be careful not to damage the mould). Inspect visually the formed element and, if necessary, correct thickness and shape. Bring the formed element using the modeling wax to the final shape. Place the formed element on casting one with wax, fasten in ring and then flood with protecting compound observing instructions of compound manufacturer. The formed element shall be cast and process in standard manner.

Danger and counteraction

Information for dental technicians:

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

P210 Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: The place of product contact with skin wash thoroughly with soap with water.

P501 Dispose of contents/container to an approved waste disposal plant.

While working with the material do not: consume meals and smoke tobacco.

Work should be carried out in good ventilated premises. Use stand ventilation.

Powder:

In case of contact of the powder particles with eyes the mechanical irritation may occur.

In case of contact acrylic powder with eyes wash it with water and contact physician.

Liquid:

Highly flammable liquid and vapor (containing methyl methacrylate), flash point +10 °C.

The place of liquid contact with skin wash thoroughly with soap water.

In case of accidental swallowing contact the physician.

Do not pour liquid into drainage system.

Storage and handling

Liquid is sensitive to influence of temperature and UV radiation under the influence of which the uncontrolled reaction of polymerization may occur.

Product components shall be stored in sales containers, in dark and good ventilated compartment, with temp. 5-25 °C, place not accessible by children.

Do not use the product after expiry date specified on the sales packaging and labels of the product ingredients.

Disposal considerations

Dispose of in accordance with applicable regulations.

Carton packaging, labels, instructions - the possibility of recycling, dispose of as communal not dangerous waste. Powder component and a packaging are not dangerous. Dispose of as communal not dangerous waste. In case of plastic container - the possibility of recycling.

Liquid and liquid packaging - hazardous waste. Pass to an authorized waste contractor.

NOTE!

Material for professional use only in technical and dental laboratories as well as dental clinics. Instructions in written form, oral or information relayed during Producers presentation result from actual level of dental technique and technology. It should be treated as obligatory information. Information's obtained from third persons (due to their different level of authorisation) do not release the user from personal check whether the preparation is intended to the designated application. Due to the fact that the application takes place without the possibility of control from the producer's side; therefore it limits the producers liability to the value of provided material.

We shall be very grateful for all comments related to products, please mail it to our address.

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