IFU APLODENT HOT	INSTRUCTIONS FOR USE – EN
-------------------------	----------------------------------

1. Product description

Aplodent Hot: Acrylic material for fast heat polymerization for dental denture bases.

Classification according to ISO 20795-1: Type 1, Class 1

Processing parameters: the data necessary for the correct processing of materials are given in [Table 1](#).

Indications for use

Aplodent Hot: Removable complete dentures (edentulous); partial dentures, indirect lining of removable dentures.

Clinical benefits

Restoration of biting and chewing functions lost as a result of partial or complete loss of dentition.

Target group(s) of patients and intended uses of devices

Users: Dental technology professionals - certified dental technicians making dental prostheses.

Patients: People of all ages and genders, with partial or complete missing teeth, eligible for prosthetic treatment with removable dentures, children, adolescents, adults.

Contraindications

Do not use the devices in patients with a diagnosed lack of resistance to allergenic agents, including the components of the device. In such cases, use a different material, suitable for the patient.

Do not use the product by pregnant women or nursing mothers in the dental technician group.

Product form

Separately: mix of powder in plastic containers and liquid in bottles.

Main components

Powder: polymethylmethacrylate (PMMA) copolymers.

Liquid: methyl methacrylate, ethylene glycol dimethacrylate.

Product compatibility

For use with acrylic artificial teeth.

Safety instructions

1 minute and leave in the closed vessel.

Warnings

H225 Highly flammable liquid and vapour.

H315 Irritating to skin.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

The components of the product and the dust generated during the manufacture of the prosthesis can cause mechanical irritation to the eyes.

Temperature- and UV-sensitive liquid. Uncontrolled polymerisation reaction possible.

Precautions

In case of any allergic reaction, stop work immediately and consult a dermatologist.

P210 Keep away from heat/sparks/open flames/hot surfaces. - Smoking prohibited.

P280 Avoid breathing dust/fume/gas/mist/vapour/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P303 + P361 + P353 IF ON SKIN (or hair); immediately remove all contaminated clothing. Rinse skin under a jet of water (or shower).

P501 Dispose of contents/container to an authorised waste collector.

It is recommended to use HDPE polyethylene gloves for hand protection, a filtering face mask suitable for liquids for respiratory protection, safety goggles and work clothing.

The use of forced room ventilation and stationary exhaust systems is advisable.

Rinse eyes with water if product component or dust comes into contact with eyes. Contact with a doctor is advisable in case of irritation.

Do not eat or drink while handling the material. In case of accidental ingestion of an ingredient of the product, contact a doctor immediately.

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Before using the product again, check that the liquid retains its properties. If there are visible signs of density change, turbidity or colour change, discontinue use of the product.

Information to be provided to patients / prosthesis users:

In order to reduce residual monomer, it is advisable to soak the denture in water during the first three days of use.

If symptoms of an allergic reaction or irritation occur, discontinue use of the denture immediately and consult a doctor.

Strong oxidising agents contained in cleaning preparations and mouthwashes can alter the colour of dentures after polymerisation. The denture should be cleaned after each meal with a soft brush under cool running water. Unscented soap, various denture cleaners and denture care products that do not contain abrasive substances can be used for cleaning.

It is advisable to wash the dentures from time to time in a vibration washers filled with a special cleaning liquid for dentures.

Wash and store the denture for the period of interruption in use in a special container designed for storing dentures.

Do not leave the denture in the sun or near sources of heat or open flame.

The prosthesis is not resistant to prolonged exposure to alcoholic beverages.

The prosthesis may change colour or lose its original physical properties if exposed to chemicals or improper use.

During use, the patient should have the prosthesis checked periodically by a dentist.

Report any damage to the denture to the dentist. Repairs may only be carried out by a dental technician.

If a piece of damaged denture is swallowed, seek medical attention immediately. The material after polymerisation is not visible on X-rays.

It is advisable to pad or rebase the dental prosthesis plate in case of atrophy of the alveolar process.

Side effects

Dental technicians

Irritation, redness or hypersensitivity symptoms including the occurrence of anaphylactic shock in case of allergy to any of the product components.

Patients

Irritation, inflammation and allergic reactions of the oral cavity caused by residual monomer leaching from the denture.

Irritation and inflammation of the mouth caused by fungi and bacteria adhering to the fabric denture.

Storage conditions

Heat- and UV-sensitive liquid. Protect from light and UV radiation.

Store the components of the product in the original sealed packaging, in a dark and well-ventilated room, at a temperature of 5-25°C and out of the reach of third parties.

Waste treatment

Dispose of in accordance with current business waste management regulations.

Cardboard packaging, labels, instructions: non-hazardous separated municipal waste - paper.

Powder: non-hazardous mixed municipal waste.

Powder packaging: non-hazardous segregated municipal waste - plastic.

Liquid and liquid packaging - hazardous waste, follow hazardous waste management regulations (storage, transport and disposal) and

hand over to a professional hazardous waste collector.

Drections for use

Preliminary steps

Read the instructions carefully before using the product.

Before making the prosthesis, ensure that the patient for whom the prosthesis is being made does not have a diagnosed lack of immunity to allergenic agents.

Each time before reusing the product, check that the fluid retains its properties.

If there are visible signs of density change, turbidity or colour change , replace the fluid with a new one.

Preparation

Teeth preparation

Use classic acrylic teeth. During setting, the teeth should have a grinded mucosal surface. Thoroughly clean the mucosal surface of the acrylic teeth from the isolation agent and moisten it with monomer. In special cases, follow the recommendations of the tooth manufacturer

Flask preparation

Use Class II gypsum (Stodent II) for the polymerization flask, **Stodent III** is particularly recommended for working models. Coat the scaled gypsum mold with a brush **twice with iso-sol** insulating liquid. Wait until the insulating film has formed.

Preparation of the material

Shake the powder and liquid containers three times, as shown on the packaging and product labels , to facilitate dispersion of the powder in the monomer.

Measuring and mixing ingredients

Recommended mixing ratio: (Table 1.). Absolutely observe the stated mixing ratio of the components. Failure to follow the instructions may result in the formation of a porous structure, a reduction in mechanical properties and an increase in the residual monomer content in the finished prosthesis.

Mixing of ingredients

Into a chemically resistant vessel with smooth inner surfaces (plastic, ceramic, porcelain, glass or stainless steel), pour the appropriate amount of liquid and carefully measure out the amount of powder used. Stir the powder and liquid vigorously with a spatula for about 1 minute and leave in the closed vessel.

Dough tying and kneading

Check from time to time whether the material is ready to be kneaded (no sticking to gloves or the sides of the bowl). The setting times for the dough are given in Table 1. The setting times refer to an ambient temperature of 23°C. Higher temperatures shorten and lower temperatures lengthen the setting time. When the material is not sticky, it should be kneaded thoroughly by kneading within 1-2 minutes. When kneading the acrylic dough, hands (gloves) must not be wet and must not be lubricated with foreign substances such as cosmetic creams, talc, etc.

Pressing

Place the prepared material in the polymerization flask and press slowly until the excess acrylic dough comes out of the flask. As an additional precaution against too little acrylic dough in the polymerisation flask, place a small additional amount of dough in the centre of the flask. Do not press the acrylic dough past its expiry time.

Polymerisation

The products belong to the group of materials that are thermally cured in an aqueous environment. They do not require the use of pressure equipment. Use standard non-pressure aqueous thermal polymerisation equipment for polymerisation.

The correct parameters for the polymerisation process are given in Table 1.

For thick dentures, use extended polymerisation over time.

Cooling and removal from the can

The cooling process parameters after polymerisation are given in Table 1.

They guarantee minimal shrinkage after polymerisation, which has an impact on the dimensional stability of the prosthesis and its retention to the palate.

Remove the frame with the flask from the polymeriser and cool to room temperature. Open the fully cooled polymerisation flask and remove the denture from the plaster.

Residual monomer reduction

Respect the specified mixing ratios, temperature and setting time of the acrylic dough and the polymerisation process parameters for the material.

In order to reduce the residual monomer content in the denture base after polymerisation, the finished denture should be soaked in water at room temperature for 12 hours before handing it over to the dentist; this step can also be performed by the patient after handing over the denture by the dentist.

Final treatment of the prosthesis

Use the standard method of treatment. Polish with **4Shine** products, **Pumice** and **Polishing Paste** by **Everall7**.

The denture plate should be polished thoroughly. The scratches left behind are the place where contaminants accumulate. Do not overheat the material during processing and polishing - permanent changes to the material's in the structure of the material leading to an increase in the residual monomer content in the finished prosthesis.

Corrections and repairs to damaged dentures

Corrections and repairs to damaged dentures can be carried out with **Aplodent Cold** acrylic material after wetting the component to be corrected or repaired with liquid.

Packaging and transportation of the denture to the dental surgery: The finished denture should be protected from mechanical damage and heat.

Hygiene, sterilisation, disinfection

Products

The devices are intended for repeated use. They are not supplied sterile nor are they required to be sterilised or disinfected prior to use. Use the principles of good dental laboratory practice to protect the devices from contamination.

Prosthesis made of materials.

There is a risk of transmission of infection. Dentures should be washed and disinfected by the dentist before being placed in the patient's mouth for the first time.

For disinfecting dentures made of the material, it is recommended to use disinfectants intended by their manufacturers for disinfecting dental prostheses, following the manufacturer's recommended procedures.

Alternatively, agents based on at least 65% ethyl or propyl alcohol can be used. The prosthesis should be sprayed with the disinfectant solution and then wiped gently dry with a paper tissue. Ready-to-use disposable wipes soaked in the disinfectant solution may also be used. In addition, the denture should be washed with warm, soapy water. Inform your dentist accordingly.

See imprint on product labels and packaging for batch number and expiry date.

The expiry date refers to the unopened packaging when stored correctly. Do not use the product after the expiry date.

Additional information

Serious incident reporting; all serious incidents involving products should be reported to the manufacturer and to the competent authority in accordance with local legislation.

Data necessary for the identification of the product and package contents

Reference number, product name, colour, pack size, batch number and expiry date - printed on the product packaging.

See table below for package contents.

Assortment *

Colour shades *

Aplodent Hot
0 - transparent
V 26 - pink veined
V 27 - milky pink veined
V 29 - coppery pink veined
T 41 - pink reddish
V 42 - pink reddish veined
V 65 - plummy pink veined

Aplodent Hot	Powder	Liquid
	4 kg, 2 kg, 1 kg, 500 g	1000 ml, 500 ml, 250 ml

* Please contact your distributor for information on product availability.

Application parameters (1)

Application parameters	Aplodent Hot	
1	Mixing ratio powder / liquid	35 g / 14 ml
2	Dough setting time (23±2°C)	approx. 10 min
3	Shelf life of the dough	approx. 25 min
4	Polymerisation	Put the flask in 100° C water, turn off the heat and leave for 15 minutes, then boil for 20 minutes.
5	Cooling	Spontaneous on air for 30´→25° C

Please send any comments on our products to our address.

For more detailed information on all Everall7 materials, please call:

+48 602 701 001

Information provided in any way, even during demonstrations, does not invalidate the instructions for use. Users are obliged to check that the product is suitable for the intended use. The manufacturer is not liable for damage, including to third parties, resulting from resulting from non-compliance with the instructions or unsuitability for use. The liability of the manufacturer is in any case limited to the value of the products supplied.

MD

This symbol indicates a medical device as defined in EU Regulation 2017/745.

IFU APLODENT HOT	MANUAL DE USO – ES
-------------------------	---------------------------

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