VILLACRYL IT	INSTRUCTIONS FOR USE – EN
---------------------	----------------------------------

patient’s mouth.

Provide your dentist with appropriate information.

Additional Notes

Keep the powder mixed with the liquid in a closed container.

The plasticity time of the acrylic dough depends on the ambient temperature. In case of any allergic reaction, stop working immediately and consult an allergist.

Application parameters (1)

Composition
Powder: polymethyl methacrylate (PMMA), ethyl methacrylate and ethyl acrylate copolymer (PMMA/EA), calcium carbonate, benzoyl peroxide (BPO), titanium dioxide, pigments

Liquid: methyl methacrylate (MMA), ethylene glycol dimethacrylate (EGDMA), N,N-dimethyl-p-toluidine (DMPT), hydroquinone.

Operating parameters: The technical data necessary for proper processing of materials are given in Table 1.

Intended use of the products

Indications: Preparing custom impression trays in the process of preparing a denture.

Contraindications

Do not use the product in persons allergic to any of the preparation components.

Acrylic monomers and benzoyl peroxide contained in the polymer are allergenic substances. In the oral environment they can be rinsed out of the material after polymerization and affect the patient’s face. Do not use the products in patients with diagnosed allergies immunodeficiency.

If you notice any allergic symptoms, stop using the device and consult with your doctor.

Do not use the product by pregnant women and nursing mothers in the group of dental technicians.

Target group(s) of patients and intended users of the devices

Users: Dental technicians making dentures, dentists taking impressions on individual trays.

Patients: Persons requiring dentures (removable complete and partial dentures, framework dentures and implant-based dentures) in situations where the use of a standard impression tray does not guarantee the required impression accuracy.

Side effects

Irritation, redness or signs of hypersensitivity, including the occurrence of anaphylactic shock, may occur in the case of allergy to any of the components of the products.

Clinical benefits

This need to reconstruct the dynamic situation of the prosthetic field in the dentures making process.

Obtaining an accurate impression of the prosthetic field in post-implantation dentures where the use of a standard impression tray does not guarantee the required impression accuracy.

Data necessary to identify the product and the contents of the package

Reference number, product name, color, package size, batch number and expiry date - imprint on the product packaging.

Package content - see table below.

Range*

Color shades

Villacryl IT
green
pink

Packaging

Villacryl IT	Powder + liquid set	Powder	Liquid
	750 g powder + 200 ml liquid	750 g	200 ml

*Please contact your distributor to obtain information about the availability of the range.

Method of Use (Step by step instructions)
Initial operations
Read the instructions carefully before using the product.
Before making the custom impression tray, make sure that the patient for whom the tray is made does not have a diagnosed allergenic immunodeficiency.
It is recommended to use gloves made of HDPE polyethylene for hand protection, a protective mask for respiratory protection, protective glasses and work clothing.

Preparing

Cost the dental plaster model twice with **izo-sol** isolating liquid. Wait until the isolator is completely dry.
Material preparation
Shake the powder and liquid containers 2 or 3 times as shown on the packaging and product labels to facilitate polymer dispersion in the liquid.
Open the powder and liquid containers as shown on the packaging and product labels.

Measuring and mixing of ingredients
Recommended mixing ratio: 21 g powder per 6 ml (5,6 g) of liquid
The stated mixing ratio must be strictly observed.

Component mixing

Pour an appropriate amount of liquid into a chemically resistant container and carefully measure the amount of powder used. Mix powder and liquid properly with spatula.
Despite the impression that the resin after the mixing is relatively dry in the first stage of preparation do not add more liquid than it is indicated in the instruction.
Making
It is recommended to use gloves from polyethylene HDPF for hands protection.
After mixing powder and liquid the acrylic dough should be immediately processed by kneading for about 1 minute.
Tray forming
After mixing the powder with the liquid, the dough remains plastic for a maximum of 4-5 minutes depending on the ambient temperature.

On the plaster model, form a tray, cut the edges to remove excess dough. Prepare the tray holder, and before connecting the holder to the tray, moisten the surface with liquid (monomer).
Material polymerization
The product belongs to the group of chemically curable materials.
The material is subject to self-polymerization at room temperature. The total polymerization time is 8 - 12 minutes from mixing the powder with the liquid.
Final treatment of the prepared impression tray
Use the standard method of treatment. The material is dusty during processing.

Packaging and transport of the impression tray to the dental office: Protect the impression tray from mechanical damage and keep it and hand it over to the dentist.

Disinfection of the tray before giving to the patient

Since there is a risk of transmission of infection, the tray should be washed and disinfected by the dentist before it is used in the patient’s mouth.
Provide your dentist with appropriate information.
Additional Notes
Keep the powder mixed with the liquid in a closed container.
The plasticity time of the acrylic dough depends on the ambient temperature. In case of any allergic reaction, stop working immediately and consult an allergist.
Application parameters (1)
Composition
Powder: polymethyl methacrylate (PMMA), ethyl methacrylate and ethyl acrylate copolymer (PMMA/EA), calcium carbonate, benzoyl peroxide (BPO), titanium dioxide, pigments
Liquid: methyl methacrylate (MMA), ethylene glycol dimethacrylate (EGDMA), N,N-dimethyl-p-toluidine (DMPT), hydroquinone.
Operating parameters: The technical data necessary for proper processing of materials are given in Table 1.
Intended use of the products
Indications: Preparing custom impression trays in the process of preparing a denture.
Contraindications
Do not use the product in persons allergic to any of the preparation components.
Acrylic monomers and benzoyl peroxide contained in the polymer are allergenic substances. In the oral environment they can be rinsed out of the material after polymerization and affect the patient’s face. Do not use the products in patients with diagnosed allergies immunodeficiency.
If you notice any allergic symptoms, stop using the device and consult with your doctor.
Do not use the product by pregnant women and nursing mothers in the group of dental technicians.
Target group(s) of patients and intended users of the devices
Users: Dental technicians making dentures, dentists taking impressions on individual trays.
Patients: Persons requiring dentures (removable complete and partial dentures, framework dentures and implant-based dentures) in situations where the use of a standard impression tray does not guarantee the required impression accuracy.
Side effects
Irritation, redness or signs of hypersensitivity, including the occurrence of anaphylactic shock, may occur in the case of allergy to any of the components of the products.
Clinical benefits
This need to reconstruct the dynamic situation of the prosthetic field in the dentures making process.
Obtaining an accurate impression of the prosthetic field in post-implantation dentures where the use of a standard impression tray does not guarantee the required impression accuracy.
Data necessary to identify the product and the contents of the package
Reference number, product name, color, package size, batch number and expiry date - imprint on the product packaging.
Package content - see table below.

Application parameters	Villacryl IT	
1	Powder / liquid mixing ratio <p>21 g powder per 6 ml (5,6 g) of liquid</p>	
2	Dough kneading time <p>1’</p>	
3	Material working time <p>4-5’</p>	
4	Polymerization <p>The material is subject to self-polymerization at room temperature.</p>	8 - 12’

Technical data (2)	
Technical Data:	Villacryl IT
Bending strength	min. 15 MPa

Warning

Information for dental technicians

Physical properties of the product: Liquid: highly flammable liquid and vapor (contains methyl methacrylate), flash point +10°C. Keep away from heat/spark/open flames/hot surfaces. - Smoking is not allowed.

Skin contact: Liquid, fine crystalline poly(methyl methacrylate powder and dusts resulting from the processing of the product may be irritating to skin. Wear protective gloves/protective clothing/eye protection/face protection.

Allergic reactions: The liquid and benzoyl peroxide contained in the polymer may cause an allergic skin reaction. IN CASE OF CONTACT WITH THE SKIN: Wash with plenty of soap and water.

Airway contact: Liquid, powder and the resulting dust can cause respiratory irritation. Avoid breathing in dust/smoke/gas/mist/vapor/sprayed liquid.

Eye contact: Accidental direct eye contact of the product’s components and dust resulting from the processing of the product may lead to eye irritation. In the event of contact with the eyes, rinse with water. Contact with an ophthalmologist is recommended in case of irritation.

Ingestion: Do not eat. In case of accidental consumption of an ingredient, contact a doctor.

Possibility of infection transmission The procedure being used enables the transmission of the infection. Apply the principles of good laboratory practice in technical and dental laboratories. Disinfect the tray before first use on the patient.

Safety measures

Information for dental technicians

Do not use the product if you are diagnosed with allergies immunodeficiency.

Standard protective measures used in technical and dental laboratories: protective masks, gloves and aprons are sufficient to safely handle the product.

In addition, it is advisable to use forced ventilation of the rooms and exhaust systems in order to reduce the risk of harmful effects of chemical substances contained in the products on the technician’s body.

In case of a contact sensitization, stop using the device and consult an allergist. In case of anaphylactic shock, take immediate rescue action.

In case of accidental consumption of an ingredient in the product, contact a specialist or poison center.

Do not use the products in patients with diagnosed immunodeficiency concerning allergens. Using the product may cause various types of sensitization reactions including anaphylactic shock. In such situations, use other material that is appropriate for the patient. Avoid inhalation of liquid vapor, direct contact of product ingredients and moist acrylic dough with skin and eyes.

Harmful to the environment. Do not pour liquids into the sewers.

Hygiene, sterilization, disinfection

The product is intended for multiple use. It is not supplied sterile and does not require sterilization or disinfection prior to use. To protect the product against contamination, apply the principles of good laboratory practices in dental laboratories.

The tray made of materials

There is a risk of transmission of infection. The tray should be cleaned and disinfected by the dentist before being placed in the patient’s mouth for the first time.

For the disinfection of trays made of materials, it is recommended to use disinfectants for disinfecting dental prostheses, following the procedures recommended by the manufacturer.

It is also possible to use agents based on at least 65% ethyl or propyl alcohol. The tray should be sprayed with a disinfectant solution, and then gently wiped dry with a paper tissue. You can also use ready-to-use disposable wipes soaked in a disinfecting solution. Additionally, the tray should be washed with warm soapy water.

Information to be given to the dentist by the technician
In case of any allergic reaction to the tray material, immediately stop taking an impression on the patient and consult a dermatologist.
If a piece of a damaged tray is swallowed while taking an impression, refer the patient to the doctor immediately. The material is not visible on X-rays.

Storing and handling the product

The expiry date refers to unopened properly stored packaging.

The product is designed for repeated use.

The liquid is sensitive to temperature and UV radiation, under the influence of which an uncontrolled polymerization reaction may occur. Each time before using the product again, check that the liquid retains its properties. If there are visible signs of a change in density, opacification or change in color, stop using the product.

Store product components (powder and liquid) in original packaging, in a dark and well-ventilated room, at a temperature of 5-25°C in a place inaccessible to third parties.

Do not use the product after the expiry date indicated on the marketing packaging and component labels.

Waste treatment

Dispose of in accordance with current waste regulations. Cardboard packaging, labels, instructions – recyclable, disposable as non-hazardous municipal waste. Powder components and powder packaging are not hazardous waste. Dispose of as non-hazardous municipal waste. In case of packaging, it is possible to be recycled. Liquid ingredients and liquid packaging – hazardous waste. Hand it over to an authorized waste collection entity.

Important points

Material to be used only in technical and dental laboratories

Information provided in any way, even during a demonstration, does not invalidate the instructions for use. Users are required to check

that the product is suitable for the intended use. The manufacturer shall not be liable for damage to third parties, resulting from failure to observe the instructions or unsuitability for use. The liability of the manufacturer is in any case limited to the value of the products supplied.

Report all serious incidents involving the medical device to the manufacturer and the relevant authorities.

Please send any comments concerning our products to our address.

MD This symbol means a medical device as defined in the EU Regulation 2017/745.

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