VILLACRYL SP	INSTRUCTIONS FOR USE - EN
---------------------	----------------------------------

I. Product description
Villacryl SP Acrylic resin for making acrylic parts in framework dentures and complete acrylic dentures with pouring technique.

ISO classification 20795-1-Type 2, Class 2.
Processing parameters: necessary parameters for the proper processing of the materials are given in **Table 1**.

Preparation: Before reusing the product, check if the fluid retains its properties. In the event of visible traces of density change, cloudiness or discoloration, stop using the product.

Indications for use
Villacryl SP: acrylic parts in framework dentures, removable acrylic dentures, repairs and relinings of dental dentures

Clinical benefits:
Restoration of biting and chewing functions lost as a result of partial or complete loss of teeth.
Target group or target groups of patients and intended uses of the devices
Users: Professionals in the dental technology industry – certified dental technicians producing dental prostheses.

Patients: People of different ages and sexes, with partial or complete lack of teeth, eligible for prosthetic treatment with removable dentures, children, adolescents, adults.
Contraindications
Do not use the product by people diagnosed with a lack of resistance to allergenic factors, including the ingredients of the product. In such situations, use a different material suitable for the patient.
Do not use the product by pregnant women and nursing mothers in the group of dental technicians.

Product form
Separate: a mixture of powder in plastic containers and liquid in bottles. The product is available in the form of a set: powder + liquid or separately as a powder and a liquid.

Chemical composition
Powder: polymethyl methacrylate (PMMA), barbituric acid derivative, pigments, viscose fibers *
Liquid: methyl methacrylate (MMA), cross-linking agent, N, N-dimethyl-p-toluidine (DMPT), hydroquinone
* Depending on the color variety, the composition of the product may slightly differ in terms of quantity and quality

Product compatibility
For use with acrylic artificial teeth and metal alloys used in dental technology industry

Safety instructions
Warnings
H225 Highly flammable liquid and vapour.
H315 Causes skin irritation.
H317 May cause an allergic skin reaction.
H335 May cause respiratory irritation.
The ingredients of the product and the dust generated during the fabrication of dentures may cause mechanical eye irritation. The liquid is sensitive to temperature and UV radiation. An uncontrolled polymerization reaction may occur.
Precautions
In case of any allergic reaction, stop work immediately and consult a dermatologist.
P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P261 Avoid breathing dust/fume/gas/mist/vapour/spray.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P302+P352 IF ON SKIN: Wash with plenty of water.
P501 Dispose of the contents/container in an authorized consignee of hand protection.
It is recommended to use gloves made of HDPE polyethylene for waste protection, a filtering fluidcompatible face mask for respiratory protection, goggles and workwear.
It is advisable to use forced ventilation of rooms and workplace exhaust systems.
In the event of contact of the product ingredient or dust with the eyes, rinse them with water. Contact a doctor is recommended in case of irritation.
Do not eat or drink while working with the material. In the event of accidental ingestion of an ingredient of the product, consult a doctor immediately.
Do not pour the liquid down the sanitation system.
Before reusing the product, check if the fluid retains its properties. In the event of visible traces of density change, cloudiness or discoloration, stop using the product.

Information to be provided to patients / prosthes users
To reduce the residual monomer, it is advisable to during the first three days, soak the prosthesis in water when not in use.
If any symptoms of an allergic reaction or irritation occur, discontinue using the prosthesis immediately and consult a doctor.
Strong oxidizing components contained in cleaning agents and mouthwashes may change the color of the prosthesis after polymerization. After each meal, the prosthesis should be washed with a soft brush under cool, running water. For cleaning, can be used unscented soap, various preparations for cleaning and caring for dental prostheses that do not contain abrasives. It is recommended to wash the dentures from time to time in a vibration washer filled with a special cleaning liquid for dentures.
When not in use, wash the denture and store it in a special denture case.
Do not leave the prosthesis in the sun, near sources of heat and open flames.
The prosthesis is not resistant to long-term effects of alcoholic beverages.
The denture may change color or lose its original physical properties when exposed to chemicals or if it is misused.
During use, the prosthesis should be periodically inspected in dental office.

Report any damage of the denture to the dentist. Repairs may be only performed by a dental technician.
If a piece of damaged dentures is swallowed, consult a doctor immediately. The material after polymerization is not visible on X-ray pictures.
It is recommended to retire or rebase the denture plate in the event of the disappearance of the alveolar process.

Side effects
Dental technicians
Irritation, redness or symptoms of hypersensitivity, including the occurrence of anaphylactic shock, in case of allergy to any of the product ingredients.
Patients
Irritation, inflammation and allergic reactions caused by residual monomer washed out of the denture.
Irritation and inflammation of the oral cavity caused by fungi and bacteria adhering to the prosthesis.

Storage conditions
Product is sensitive to temperature and UV radiation. Protect against light and UV radiation.
Store the product components in their original, tightly closed containers, in a dark and well-ventilated room, at a temperature of 5-25 °C, out of the reach of children.
Waste disposal
Dispose in accordance to applicable regulations of economic activity waste management.
Cardboard packages, labels, instructions - non-hazardous segregated municipal waste – paper.
Powder: non-hazardous mixed municipal waste.
Powder packaging: non-hazardous segregated municipal waste – plastic
The liquid and its packaging - hazardous waste, comply with the regulations on hazardous waste management (storage, transport and disposal) and hand it over to a professional hazardous waste consignee.

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Additional information
Reporting of serious incidents: Any serious incidents related to the products must be reported to the manufacturer and the competent authority in accordance with local regulations.
Data necessary to identify the product and the contents of the package
Reference number, product name, color, package size, batch number and expiry date - imprint on the product packaging.
Package content - see table below.

Assortment*
Color shades
Villacryl SP
0 - transparent
V2 – milk pink veined
V4 – pink veined

	Powder + liquid kit	Powder	Liquid
Villacryl SP	500g powder + 300 ml liquid;	500 g	300 ml

* Information on the availability of the assortment is provided by the distributor.
Processing parameters (1)

Application parameters	Frameworks	Pouring technique in agar
1 Mixing ratio powder/liquid	10 g / 6,7 g (7 ml)	10 g / 5 g (5,2 ml)
2 Mixing time (23±2°C)	30” - 60”	30” - 60”
3 Polymerization	min. 20” 50°C - 60°C Pressure 2 bar	30” 60°C Pressure 2 bar
4 Cooling	Remove the model from the polymerizer and cool to room temperature.	Remove the flask from the polymerizer and cool to room temperature.

Problem prevention (2)
See website: **www.everall.pl**

Please send any comments regarding our products to our address.

Information provided in any way, even during the demonstration, does not invalidate the instructions for use. Users are responsible for checking that the product is suitable for the intended use. The manufacturer is not liable for damages, including to third parties, resulting from non-compliance with the instructions or non-use. The manufacturer’s liability is in any case limited to the value of the products delivered.

MD	The symbol means a medical device as defined in EU Regulation 2017/745.
-----------	---

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Informacje dodatkowe
Reporting of serious incidents: Any serious incidents related to the products must be reported to the manufacturer and the competent authority in accordance with local regulations.
Data necessary to identify the product and the contents of the package
Reference number, product name, color, package size, batch number and expiry date - imprint on the product packaging.
Package content - see table below.

Assortment*
Color shades
Villacryl SP
0 - transparent
V2 – milk pink veined
V4 – pink veined

	Powder + liquid kit	Powder	Liquid
Villacryl SP	500g powder + 300 ml liquid;	500 g	300 ml

* Information on the availability of the assortment is provided by the distributor.
Processing parameters (1)

Application parameters	Frameworks	Pouring technique in agar
1 Mixing ratio powder/liquid	10 g / 6,7 g (7 ml)	10 g / 5 g (5,2 ml)
2 Mixing time (23±2°C)	30” - 60”	30” - 60”
3 Polymerization	min. 20” 50°C - 60°C Pressure 2 bar	30” 60°C Pressure 2 bar
4 Cooling	Remove the model from the polymerizer and cool to room temperature.	Remove the flask from the polymerizer and cool to room temperature.

Problem prevention (2)
See website: **www.everall.pl**

Please send any comments regarding our products to our address.

Information provided in any way, even during the demonstration, does not invalidate the instructions for use. Users are responsible for checking that the product is suitable for the intended use. The manufacturer is not liable for damages, including to third parties, resulting from non-compliance with the instructions or non-use. The manufacturer’s liability is in any case limited to the value of the products delivered.

MD	The symbol means a medical device as defined in EU Regulation 2017/745.
-----------	---

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