|               |                          |
|---------------|--------------------------|
| VILLACRYL STC | INSTRUCTION FOR USE - EN |
|---------------|--------------------------|

#### Product Characteristics

Self-curing acrylic resin for temporary crowns and bridges and for repairing temporary acrylic-faced crowns and bridges. Material in powder and liquid form.

Preforming of the crown and curing of the resin takes place directly in the patient's oral cavity.

#### Composition

#### Villacryl STC

Powder: polymethyl methacrylate (PEMA), 1-benzyl-5-phenyl-barbituric acid (BPBA), benzoyl peroxide (BPO), titanium dioxide, pigments.

Liquid: methyl methacrylate (MMA), trietylene glycol dimethacrylate (TEGDMA), hexanediol 1,6-dimethacrylate (1,6 HDMA), N,N-dimethyl-p-toluidine (DMPT), hydroquinone (HQ).

**Operating parameters:** The technical data necessary for the proper processing of materials are given in the text and **Table 1**.

#### Intended use of the product

#### Indications

Making temporary crowns and bridges for a transitional period until the permanent filling is completed.

#### Contraindications:

Do not use the products in patients with diagnosed immunodeficiency concerning allergens.

Do not use the material to make / repair temporary crowns for more than 30 days.

Do not use the device in clinical applications other than those specified in the instructions.

Do not use the product by pregnant women and nursing mothers in the group of dental technicians.

#### Target group(s) of patients and intended users of the devices

**Users:** Dentists making temporary crowns

**Patients:** Persons of various ages, with small dentition losses allowing the use of acrylic crowns or bridges, during the period of prosthetic treatment until the permanent restoration is made.

#### Side effects

During crown formation in the oral cavity, the patient may suffer from the effects of monomer fumes on the mucous membrane and soft tissues of the oral cavity, as a slight burning sensation or unpleasant taste or smell. This can lead to vomiting reflexes in sensitive patients. These sensations decrease with time and the binding of the material until they disappear completely.

Irritation, redness or signs of hypersensitivity, including the occurrence of anaphylactic shock, may occur in the case of allergy to any of the components of the products.

Oral irritation and inflammation caused by fungi and bacteria adhering to a crown or bridge made of material is possible in the absence of proper hygiene.

#### Clinical benefits

Restoration of the lost biting and chewing functions as a result of partial loss of dentition until the time of permanent denture restoration.

Improving the appearance, comfort and well-being during the treatment period and waiting for the permanent restoration of people who have lost their natural teeth.

Repair of temporary acrylic veneer crowns damaged during use.

Improving the comfort and well-being of users of acrylic crowns that have been mechanically damaged during use.

#### Data necessary to identify the product and the contents of the package

Reference number, product name, color, package size, batch number and expiry date - imprint on the product packaging.

Package content - see table below.

#### Range\*

#### Color shades

|                      |   |   |   |
|----------------------|---|---|---|
| <b>Villacryl STC</b> | 1 | 2 | 3 |
|----------------------|---|---|---|

|                      |   |               |               |                 |  |
|----------------------|---|---------------|---------------|-----------------|--|
| <b>Packaging</b>     |   |               |               |                 |  |
|                      | <b>Powder + liquid set</b>  | <b>Powder</b> | <b>Liquid</b> | <b>Hardener</b> |  |
| <b>Villacryl STC</b> | <b>Villacryl STC powder 3 x 20 g + liquid 40 ml + hardener 40 g</b> | 20 g          | 40 ml         | 40 g            |  |

\*Please contact your distributor to obtain information about the availability of the range.

#### Method of Use (Step by step instructions)

#### Initial operations:

Read the instructions carefully before using the product for the first time.

The formation of the crown or bridge and its preliminary curing takes place in the patient's oral cavity. Before using the material, make sure that the patient is not diagnosed with an immunodeficiency concerning allergens or is not sensitive to taste or smell. Otherwise, the use of a different material is recommended.

Wear gloves, a protective mask and protective clothing that comply with dental practice standards and procedures to protect the health and safety of medical personnel and patients.

#### Preparing

#### Tooth preparation

Take a dental impression of the patient's oral cavity that reconstructs the condition before grinding the teeth, or prepare a standard mold suitable for the size of the tooth.

Prepare the teeth for the crown by grinding them according to the current level of knowledge and technology for the treatment process.

#### Material preparation

Select a suitable material color similar to the natural color of patient's teeth.

Shake the powder and liquid containers 2 or 3 times to facilitate powder dispersion in the monomer.

#### Measuring and mixing of ingredients

#### Recommended mixing ratio:

The stated mixing ratio must be strictly observed. Failure to comply with the stated values may lead to the formation of a porous crown structure, a reduction in its mechanical properties and an increase in the residual methyl methacrylate monomer content in the finished restoration.

#### Component mixing

Mix exactly 2 g of powder with 1 g (~1 mL) of liquid in a chemically resistant container and pour into a prepared standard mold or impression taken before grinding the tooth. The mixing time should not exceed 1 minute.

#### Crown forming

When forming the crown, insulate the ground teeth with Vaseline.

Place a mold or impression tray with the material in the patient's mouth for 4 minutes. During this time, the patient may experience slight discomfort from the effects of monomer fumes on the mucous membrane and soft tissues of the oral cavity, as a slight burning sensation or unpleasant taste or smell sensations. These sensations decrease with time and the binding of the material until they disappear completely. Avoid direct contact of moist acrylic dough with the tongue and mucous membrane as far as possible. After curing, remove the crown and remove any excess material. Instruct the patient to rinse his mouth thoroughly with water.

#### Polymerization

Place the formed and pre-cured crown in a container with hot water at 50-65°C for 15 minutes, adding the 0.5 scoop (2 g) of the Hardener component per glass of water (200 ml).

The parameters of the polymerization process guarantee a crown with appropriate parameters and ensure high patient safety associated with its use.

#### Final treatment of the crown

Process the crown and polish it thoroughly with e.g. Everall7 polishing paste. Scratches left in the denture are a place of dirt accumulation.

Do not overheat the material during the processing and polishing process – permanent changes in the material structure may occur, leading to an increase in the presence of a residual monomer in the finished crown.

Place the crown in the patient's oral cavity and place it on the ground teeth using cement for temporary crowns.

#### Corrections and repairs of damaged crown

Corrections and repairs of a damaged temporary crown or bridge should be performed with Villacryl STC mass using standard methods.

#### Hygiene, sterilization, disinfection

#### Product

The product is intended for multiple use. It is not supplied sterile and does not require sterilization or disinfection prior to use. To protect the product against contamination, apply the principles of good laboratory practices in dental laboratories.

#### The crowns made of materials

There is a risk of transmission of infection. The crown should be cleaned and disinfected by the dentist before being placed in the patient's mouth for the first time.

For the disinfection of crown made of materials, it is recommended to use disinfectants for disinfecting dental prostheses, following the procedures recommended by the manufacturer.

It is also possible to use agents based on at least 65% ethyl or propyl alcohol. The crown should be grayed with a disinfectant solution, and then gently wiped dry with a paper tissue. You can also use ready-to-use disposable wipes soaked in a disinfecting solution. Additionally, the crown should be washed with warm soapy water. Provide the appropriate information to the dentist.

#### Information to be provided to patients / prosthesis users

If any symptoms of an allergic reaction or irritation occur, discontinue using the crown immediately and consult a doctor.

Strong oxidizing compounds contained in cleaning agents and mouthwashes may change the color of the material after polymerization.

For the hygiene of crowns and bridges, popular and generally available toothpastes, threads, interdental brushes can be used, it is recommended to use a pressure device of the WATERPIK type.

The material is not resistant to long-term effects of alcoholic beverages.

The material may change color or lose its original physical properties when exposed to chemicals if it is misused.

Report any damage of the material to the dentist. Repairs may be only performed by a dental technician.

If a piece of damaged crown is swallowed, consult a doctor immediately.

The material after polymerization is not visible on X-ray pictures.

#### Additional Notes

In case of any allergic reaction, stop working immediately and consult your dermatologist.

#### APPLICATION PARAMETERS

|  |                      |
|--|----------------------|
| <b>Application parameters</b>                          | <b>Villacryl STC</b> |
| 1 Powder / liquid mixing ratio                         | 2 g / 1 ml           |
| 2 Dough preparation time (23±2°C)                      | 30 s                 |
| 3 Shelf life of the dough                              | 1 min 30 s           |
| 4 Curing time in oral cavity                           | 4 min                |
| 5 Maximum temperature during curing in the oral cavity | 37°C                 |
| 6 Polymerization                                       | 15 min, 50-65°C      |

#### Warning

**Skin contact:** The liquid, fine crystalline polymethyl methacrylate powder and dusts resulting from the processing of the product can irritate the skin. Wear protective gloves/protective clothing/eye protection/face protection.

**Allergic reactions:** The liquid, benzoyl peroxide contained in the polymer, and the dust produced during processing and grinding of the crown may cause an allergic skin reaction. IN CASE OF CONTACT WITH THE SKIN: Wash with plenty of soap and water.

**Airway contact:** liquid, powder and the resulting dust can cause respiratory irritation. Avoid breathing in dust/smoke/gas/mist/vapor/sprayed liquid.

**Eye contact:** Accidental direct contact of the product components and the dust generated by the crown with the eyes may lead to their irritation. In the event of contact with the eyes, rinse with water. Contact with an ophthalmologist is recommended in case of irritation.

**Ingestion:** Do not eat. In case of accidental consumption of an ingredient, contact a doctor.

**Possibility of infection transmission:** The procedure being used enables the transmission of the infection. Apply good laboratory practice in dental offices when making or repairing a crown or bridge. The crown or bridge should be washed and disinfected before being placed in the patient's oral cavity.

#### Safety measures

#### Dentist:

Do not eat the product if you are diagnosed with immunodeficiency concerning allergens.

Do not eat or smoke when working with the material.

Standard protective equipment: protective masks, gloves and aprons are sufficient for safe handling of the product.

In case of a contact sensitization, stop using the device and consult an allergist. In case of anaphylactic shock, take immediate rescue action.

In case of accidental consumption of an ingredient in the product, contact a specialist or poison center.

Do not pour liquids into the sewers.

**Patient:**

Patients with diagnosed immunodeficiency concerning allergens: using the product may cause various types of sensitization reactions including anaphylactic shock. In such situations, use other material that is appropriate for the patient.

Patients sensitive to taste and smell – making a crown in the oral cavity can lead to vomiting reflexes. If the reaction is strong, stop making the crown and use other material suitable for the patient.

#### Other patients:

Insulate the ground teeth with Vaseline.

When forming the crown, avoid direct contact of excess moist acrylic dough with the tongue and mucous membrane of the patient as far as possible. In case of unintentional contact of excess moist acrylic dough with the patient's mucous membrane or tongue, remove it from the oral cavity and then ask to rinse the mouth with water. If irritation or allergy symptoms occur, it is necessary to refer the patient to a specialist doctor for appropriate treatment.

In case of accidental ingestion of the material by the patient during crown forming, refer the patient to a specialist doctor or contact the poison center.

After the crown has been formed, ask the patient to rinse his mouth thoroughly with water.

**Reporting of serious incidents:** Any serious incidents related to the products must be reported to the manufacturer and the competent authority in accordance with local regulations.

#### Storing and handling the product

The expiry date refers to unopened properly stored packaging.

The product is designed for repeated use.

The liquid is sensitive to temperature and UV radiation, under the influence of which an uncontrolled polymerization reaction may occur. Each time before using the product again, check that the liquid retains its properties. If there are visible signs of a change in density, opacification or change in color, stop using the product.

Store product components (powder and liquid) in original packaging, in a dark and well-ventilated room, at a temperature of 5-25°C in a place inaccessible to third parties.

Do not use the product after the expiry date stated on the marketing packaging and labels.

#### Waste treatment

Dispose in accordance with applicable regulations of economic activity waste management. Cardboard packages, labels, instructions: non-hazardous segregated municipal waste- paper. Powder: non-hazardous mixed municipal waste.

Powder packaging: non-hazardous segregated municipal waste - plastic.

The liquid and its packaging - hazardous waste, comply with the regulations on hazardous waste management (storage, transport and disposal) and hand it over to a professional hazardous waste consignee.

After removing from the patient's oral cavity, the worn crown shell be treated as medical waste. Dispose of according to waste regulations.

#### Problem prevention

See website: **www.everall7.pl**

#### Important points

Information provided in any way, even during a demonstration, does not invalidate the instructions for use. Users are required to check that the product is suitable for the intended use. The manufacturer shall not be liable for damage, including damage to third parties, resulting from failure to observe the instructions or unsuitability for use. The liability of the manufacturer is in any case limited to the value of the products supplied.

#### MD

Please send any comments concerning our products to our address.

#### MD

This symbol means a medical device as defined in the EU Regulation 2017/745.

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