<b>IFU APLODENT COLD</b>	<b>INSTRUCTIONS FOR USE – EN</b>
--------------------------	----------------------------------

Each time before reusing the product, check that the liquid retains its properties. Stop using the product if there are visible signs of density change, turbidity or colour change.

#### 1. Product description

**Aplodent Cold:** acrylic resin for fast cold polymerisation for the repair of dental dentures.

**Classification according to ISO 20795-1:** Tipo 2, Clase 1

**Performance parameters:** the data necessary for the correct processing of the materials are given in Table 1.

#### Indications for use

**Aplodent Cold:** repair of removable complete dentures, partial dentures, acrylic frame denture parts and acrylic orthodontic appliances.

#### Clinical benefits

Restoring biting and chewing function to dentures damaged during use.

Restoration of orthodontic appliances damaged during use to correct the dentition or prevent secondary movement of the corrected teeth.

#### Target group(s) of patients and potential users of devices

**Users:** Dental technology professionals - certified dental technicians making dental prostheses.

**Patients:** People of all ages and genders, dental prosthesis users.

Acrylic braces users, children, adolescents, adults.

#### Contraindications

Do not use the device in patients with a diagnosed lack of resistance to allergenic agents, including the components of the device. In such cases, use a different material, suitable for the patient.

Do not carry out repairs, linings or other direct work in the patient's mouth.

Do not use the product by pregnant women or nursing mothers in the dental technician group.

#### Product form

Separately: mix of powder in plastic containers and liquid in glass bottles.

#### Main components

Powder: poly(methylmethacrylate).

Liquid: methyl methacrylate, ethylene glycol dimethacrylate, N,N dimethyl-p-toluidine.

#### Product compatibility

For use with acrylic artificial teeth and acrylic dental prostheses.

#### Safety instructions

#### Warnings

H225 Highly flammable liquid and vapour.

H315 Irritating to skin.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

The components of the product and the dust generated during denture repair can cause mechanical irritation to the eyes.

Temperature- and UV-sensitive liquid. Uncontrolled polymerisation reaction possible.

#### Precautions

In case of any allergic reaction, stop work immediately and consult a dermatologist.

P210 Keep away from heat/sparks/open flames/hot surfaces. - Smoking prohibited.

P261 Avoid breathing dust/fume/gas/mist/vapour/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P303 + P361 + P353 - IF ON SKIN (or hair): immediately remove all contaminated clothing. Rinse skin under a jet of water (or shower).

P501 Dispose of contents/container to an authorised waste collector.

It is recommended to use HDPE polyethylene gloves for hand protection; a filtering face mask suitable for liquids for respiratory protection, safety goggles and work clothing.

The use of forced room ventilation and stationary exhaust systems is advisable.

Rinse eyes with water if product component or dust comes into contact with eyes. Contact with a doctor is advisable in case of irritation.

Do not eat or drink while handling the material. In case of accidental ingestion of an ingredient of the product, contact a doctor immediately.

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Before using the product again, check that the liquid retains its properties. Stop using the product if there are visible signs of density change, turbidity or colour change.

#### Information to be given to patients / prosthesis users

In order to reduce residual monomer, it is advisable to soak the denture in water during the first three days of use.

If symptoms of an allergic reaction or irritation occur, discontinue use of the denture immediately and consult a doctor.

The strong oxidising agents in cleaning preparations and mouthwashes can change the colour of the material after polymerisation.

Unscented soap, various denture cleaners and denture care products that do not contain abrasive substances can be used to clean the repaired denture material.

It is advisable to wash the dentures from time to time in a vibration wash filled with a special cleaning liquid for dentures.

The material is not resistant to prolonged exposure to alcoholic beverages.

The site of the repaired prosthesis may change colour or lose its original physical properties if it comes into contact with chemicals or is not properly handled.

Report any repeated damage to the denture during the period of use to the dentist. Repairs may only be carried out by a dental technician.

If a piece of damaged repaired denture is swallowed, see a doctor immediately. The material after polymerisation is not visible on X-rays.

#### Side effects

#### Dental technicians

Irritation, redness or hypersensitivity symptoms including the occurrence of anaphylactic shock. In case of allergy to any of the product components.

#### Patients

Irritation, inflammation and allergic reactions of the oral cavity caused by residual monomer leaching from the repaired denture.

#### Storage conditions

Heat- and UV-sensitive liquid. Protect from light and UV radiation.

Store the product in original, tightly closed packaging, in a dark and well-ventilated room, at a temperature of 5-25°C, out of the reach of third parties.

#### Waste treatment

Dispose of in accordance with current business waste management regulations.

Cardboard packaging, labels, instructions: non-hazardous segregated municipal waste - paper.

Powder: non-hazardous - mixed municipal waste.

Powder packaging: non-hazardous segregated municipal waste - plastic.

Liquid and liquid bottle - hazardous waste, follow hazardous waste management regulations (storage, transport and disposal) and hand-over to a professional hazardous waste collector.

#### Directions for use

#### Preliminary steps

Read the instructions carefully before using the product.

Ensure that the patient for whom the denture repair is being carried out has not been diagnosed with an immune deficiency to allergic agents.

+48 602 701 001

Information provided in any way, even during demonstrations, does not invalidate the instructions for use. Users are obliged to check that the product is suitable for the intended use. The manufacturer is not liable for damage, including to third parties, resulting from

resulting from non-compliance with the instructions or unsuitability for use. The liability of the manufacturer is in any case limited to the value of the products supplied.

<b>MD</b>	This symbol indicates a medical device as defined in EU Regulation 2017/745.
-----------	--

<b>IFU APLODENT COLD</b>	<b>MANUAL DE USO – ES</b>
--------------------------	---------------------------

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