Universal metal masking material

The product characteristics: Acrylic resin for masking metal elements in dentures, crowns, and bridges.

Classification in accordance with EN ISO 10477: Type 2, Class 1

Indications for use: skeletal dentures, crowns, and bridges faced with acrylic.

Clinical benefits

Restoration of biting and chewing functions lost as a result of partial or complete loss of teeth.

Target group or target groups of patients and intended users of the devices

Users: Professionals in the dental technology industry - certified dental technicians producing dentures, crowns and bridges.

Patients: People of different ages and sexes, with partial or complete lack of teeth, eligible for use of dentures, acrylic or bridges.

Contraindications: Do not use this product on persons allergic to any of thecomponents of the product. Do not use the product by pregnant women and nursing mothers in the group of dental technicians.

Side effects

Dental technicians

Irritation, redness or symptoms of hypersensitivity, including the occurrence of anaphylactic shock, in case of allergy to any of the product ingredients.

Patients

No side effect are known at the moment.

Data necessary to identify the product and the contents of the package

Reference number, product name, color, package size, batch number and expiry date - imprint on the product packaging.

Package content - see Table below.

REF	Set
V210Z10	Villacryl Opaker 4x7 g powder (grey, yellow, brown, pink colours) + 2x12 ml liquid + 12 ml Primer
V210Z12	Villacryl Opaker 2x7 g powder (pink colour) + 12 ml liquid + 12 ml Primer
V210L03	Villacryl Opaker Liquid 12 ml
V210L04	Villacryl Opaker Primer 12 ml

Chemical composition

Powder: acrylic polymer, benzoyl peroxide, pigments.

Liquid: acrylic monomers, N,N dimethyl-p-toluidine.

Primer: silane, isopropyl alcohol.

Product compatibility

For use with dental metals and alloys as well as acrylic materials for the preparing of dentures, crowns and dental bridges.

Method of use

Preparation: Cleanse the metal surface of any leftovers of the protective mass or metal oxides using a prosthetic sandblaster. Prepare the cleansed surface by sandblasting it with pure aluminium oxide with the granulation of 250 µm and under the pressure of 4 bar. Clean the prepared surface with acetone, apply **Opaker Primer**, and spread it evenly with compressed air (do not use oil compressors). **Opaker Primer** is used for the chemical joining of acryl-metal, and it prevents the creation of the so-called "randkluft."

Mixing: Pour a few drops of the **Opaker Liquid** monomer to a separate container, and add some **Opaker Powder** to saturate it. Stir for about 30 seconds till you get a homogenous consistency. The **Opaker Powder** powders can be mixed with one another in order to use the required colour.

Applying: Spread on the surface of the metal with a thin, opaque layer of the prepared product with the help of a small brush, and wait for about 5 minutes.

Curing: Curing the opaque product can be also performed in a stream of hot air over flame of a blowpipe. One ought be careful not to burn the product layer. Otherwise, it will become chalky, and it will not adhere to the metal surface.

If the product isured in the right way, it keeps a permanent in colour; it adheres to metal; and it covers it completely.

Treatment and finishing: After curing, wash the Opaker surface with alcohol in order to remove the unpolymerized layer. Material cured in water has a brighter color than cured in a stream of hot air.

Warnings

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes eye irritation.

H335 May cause irritation of the respiratory tract.

H336 May cause drowsiness or dizziness.

Precautions

P210 Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

P271 Avoid breathing in dust/smoke/gas/mist/vapours/spray.

P271 Use only outdoors or in a well-ventilated area.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: Wash thoroughly with soap with water.

P304+P340 IF INHALED: Take the victim to fresh air and ensure conditions for rest in the position allowing for free breathing.

P305+P351+P338 IF IN EYES: In case of contact with the acrylic powder or liquid with eyes, wash them with water and contact a doctor.
Remove contact lenses if present and easy to do - continue rinsing.
P501 Dispose of the contents/container to an approved waste disposal plant. An inflammable Primer (contains isopropyl alcohol) - flash point +15 °C. An inflammable liquid (contains) - flash point +10 °C. Keep away from fire and radiation sources.
While working with the matenal, do not consume meals or smoke tobacco.
Avoid breathing the vapour of the liquid, and the direct contact of the product components or the wet acrylic dough with the skin and eyes.
Work should be carried out in well ventilated premises.
Use workplace ventilation.
Wash thoroughly the areas of contact of the product components with skin with soap and water. In case of accidental consumption, contact the physician. Do not pour into the drainage system.

Hygiene, sterilization, disinfection

Product

The product is intended for multiple use. It is not supplied sterile and does not require sterilization or disinfection prior to use. To protect the product against contamination, apply the principles of good laboratory practices in dental laboratories.

Denture framework / metal element of the crown masked by material

It does not require sterilization or disinfection before applying the correct acrylic layer of the restoration.

Information to be provided to patients / prosthesis users:

Not applicable

Storage and handling

Product components are to be stored in original containers, in a dark and well ventilated room with temp. 5-25 °C, in a place not accessible to children. Do not use the product after the expiry date. Unused material and packaging are dangerous for the environment. Remove according to local regulations.

Disposal considerations

Dispose of in accordance with applicable regulations. Carton packaging, labels, instructions - the possibility of recycling; dispose of as communal waste not constituting a danger. The powder component and powder packaging are not dangerous waste.

Dispose of as communal, not dangerous waste. In case of packaging- the possibility of recycling. The liquid and liquid packaging - hazardous waste. Pass on to an authorized waste contractor. The primer and primer packaging - hazardous waste. Pass on to an authorized waste contractor.

Problem prevention (2)

See website: **www.everall7.pl**

Additional information

Reporting of serious incidents: Any serious incidents related to the products must be reported to the manufacturer and the competent authority in accordance with local regulations.

Instructions in the written and oral form, and information relayed during a presentation by the Producer result from the current level of dental technique and technology. They ought to be treated as applicable information. The information obtained from third persons (due to their different level of aforisation) does not release the user from checking personally whether the product is intended for the designated application. Due to the fact that the application takes place without the possibility of control on the producer’s part, it limits the producer’s liability to the value of provided material.

<p>We kindly request all the comments connected with our products to be sent to our address.</p>
<p>VILLACRYL OPAKER</p>
<p>INSTRUKCJA UŻYCIA - PL</p>

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