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**I. Product description**

**Villocryl Ortho:** Acrylic resin for making removable orthodontic appliances.  
**Villocryl Ortho Mix:** Acrylic resin for making removable orthodontic appliances.

ISO classification 2079-2: Type 1

**Processing parameters:** necessary parameters for the proper processing of the materials are given in **Table 1**.

**Indications for use**

**Villocryl Ortho:** Removable orthodontic appliances, repairs of orthodontic appliances.

**Villocryl Ortho Mix:** Removable orthodontic appliances, repairs of orthodontic appliances.

**Clinical Benefits**

Correcting minor changes in children, adolescents and adult's dentition.

Maintaining corrected dentition with orthodontic appliances. Preventive of returning teeth to the state before treatment.

Restoring the proper functions to orthodontic appliances damaged during use.

**Target group or target groups of patients and intended users of the devices**

**Users:** Professionals in the dental technology industry - certified dental technicians producing removable orthodontic appliances.

**Patients:** People of different ages and sexes, with predominance of children and adolescents, with slight malocclusion, not eligible for orthodontic treatment with fixed appliances.

People of various ages and sexes, with a predominance of children and adolescents after completed orthodontic treatment with fixed or removable appliances.

**Contraindications**

Do not use the products by people diagnosed with a lack of resistance to allergenic factors, including the ingredients of the product. In such situations, use a different material suitable for the patient.

Do not use the product by pregnant women and nursing mothers in the group of dental technicians.

**Product form**

Separate: a mixture of powder in plastic containers and liquid, color concentrates in bottles.

The product is available in the form of a set: powder + liquid + color concentrates or separately as a powder + liquid and a color concentrate.

**Chemical composition**

**Villocryl Ortho**

Powder: polyethyl methacrylate (PEMA), barbituric acid derivative

Fluid: methyl methacrylate (MMA), cross-linking agent, stabilizer, amine catalyst

Color concentrates: methyl methacrylate (MMA), oil-soluble pigments

**Villocryl Ortho Mix**

Powder: polymethyl methacrylate (PMMA)/polyethyl methacrylate (PEMA), barbituric acid derivative.

Fluid: methyl methacrylate (MMA), cross-linking agent, stabilizer, amine catalyst

Product compatibility

For use with standard orthodontic wires and screws

**Safety instructions**

**Warnings**

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

The ingredients of the product and the dust generated during the fabrication of dentures may cause mechanical eye irritation.

The fluid is sensitive to temperature and UV radiation. An uncontrolled polymerization process may occur.

**Precautions**

In case of any allergic reaction, stop work immediately and consult a dermatologist.

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: Wash with plenty of water.

P501 Dispose of the contents/container to an authorized consignee of the waste.

It is recommended to use gloves made of HDPE polyethylene for hand protection, a filtering fluidcompatible face mask for respiratory protection, goggles and workwear.

It is advisable to use forced ventilation of rooms and workplace exhaust systems.

In the event of contact of the product ingredient or dust with the eyes, rinse them with water. Contact a doctor is recommended in case of irritation.

Do not eat or drink while working with the material. In the event of accidental ingestion of an ingredient of the product, consult a doctor immediately.

Do not pour the liquid down the sanitation system.

Before reusing the product, check if the fluid retains its properties. In the event of visible traces of density change, cloudiness or discoloration, stop using the product.

**Information to be provided to patients / orthodontic appliances users**

It is advisable during the first three days, to soak the orthodontic appliance in water when not in use.

If any symptoms of an allergic reaction or irritation occur, discontinue using the orthodontic appliance immediately and consult a doctor.

Strong oxidizing compounds contained in cleaning agents and mouthwashes may change the color of the orthodontic appliance after polymerization. After each meal, the prosthesis should be washed with a soft brush under cool, running water. For cleaning, can be used unscented soap, various preparations for cleaning and caring for orthodontic appliance that do not contain abrasives. It is recommended to wash the orthodontic appliance from time to time in a vibration wash filtered with a special cleaning liquid for dentures.

When not in use, wash the orthodontic appliance and store it in a special denture case.

Do not leave the orthodontic appliance in the sun, near sources of heat and open flames.

The orthodontic appliance may change color or lose its original physical properties when exposed to chemicals or if it is misused.

During use, the orthodontic appliance should be periodically inspected in dental office.

Report any damage of the denture to the dentist. Repair may be only performed by a dental technician.

If a piece of damaged orthodontic appliance is swallowed, consult a doctor immediately. The material after application is not visible on X-ray pictures.

**Side effects**

**Dental technicians**

Irritation, redness or symptoms of hypersensitivity, including the occurrence of anaphylactic shock, in case of allergy to any of the product ingredients.

**Patients**

Irritation, inflammation and allergic reactions caused by residual monomer washed out of the orthodontic appliance.

Additionaly, the inflammation of the oral cavity caused by fungi and bacteria adhering to the orthodontic appliance.

**Storage conditions**

Product and color concentrates are sensitive to temperature and UV radiation. Protect against light and UV radiation.

Store the product components in their original, tightly closed containers, in a dark and well-ventilated room, at a temperature of 5-25 ° C out of the reach of third parties.

**Waste disposal**

Dispose in accordance with applicable regulations of economic activity waste management.

Cardboard packages, labels, instructions: non-hazardous segregated municipal waste - paper.

Powder: non-hazardous mixed municipal waste.

Powder packaging: non-hazardous segregated municipal waste - plastic.

The liquid and packaging of liquid and color concentrates - hazardous waste, comply with the regulations on hazardous waste management (storage, transport and disposal) and hand it over to a professional hazardous waste consignee.

**How to use (Step by step instructions)****Initial steps**

Read the instructions carefully before using the product.

Before making the orthodontic appliance, make sure the patient for whom the prosthesis is made has not been diagnosed with the lack of resistance to allergenic factors.

Before reusing the product, check if the fluid retains its properties. In the event of visible traces of density change, cloudiness or discoloration, stop using the product.

**Preparation**

Shake the powder and liquid containers 3 times as shown on the package and product labels to help disperse the powder in the monomer.

Open the powder and liquid containers as shown on the packaging and product labels.

**Model preparation:** glue with wax or other material (if **Form Plast**) the previously bent and prepared metal elements (clasp, bows) on the vestibular side of the dried stone model. Attach the orthodontic screw. Then use a brush to coat the plaster model **twice** with **Isol-soil** - isolating liquid. Wait for the isolation membrane to form. Soak the model base in water until completely saturated (with water).

For working models, use 3rd and 4th class stone, especially **Stodent III**, **Stodent III Ortho** and **Stodent IV**.

**Application of the fluid solution for work**

**Recommended proportion to obtain the color:** 2-4 drops of the desired color concentrate per 1 ml of liquid.

More intense color can be obtained by adding more concentrate. Mix the liquid and concentrate in an appropriate container. The solution is ready for immediate use.

**Polymerization:** Use pressure polymerization in polymerization pot filled with water. The correct parameters of the polymerization process for the material are given in Table 1.

The product is not intended for use with techniques such as, for example, curing with light, microwave or injection method.

**I. Making orthodontic appliance with salt and pepper technique (Villocryl Ortho / Villocryl Ortho Mix)**

Application onto the model: pour a thin layer of polymer on the surface of the isolated stone model. Using a pipette, moisten the powder with the liquid containing the color concentrate. Apply another portion of the powder on the damp surface and moisten it with the liquid. Repeat this operation several times until the desired thickness of the plate of the orthodontic appliance is obtained.

Polymerization: Place the stone model with formed orthodontic appliance in a pressure polymerizer.

Polymerize in water maintaining the parameters given in Table 1.

**II. Making orthodontic appliance with acrylic dough technique (Villocryl Ortho Mix)**

**Measuring and mixing**

**Recommended mixing ratio:** (Table 1). It is essential to follow the stated mixing ratio of the components.

Not following the instructions may result in the formation of a porous structure, a reduction in mechanical properties and increase the level of residual monomer.

**Mixing:** Prepare the ingredients of the product in the amount necessary to make the orthodontic appliance, minding the given mixing ratio. Pour the appropriate amount of colored solution liquid into a chemically resistant vessel and carefully measured amount of powder.

Mix powder and liquid thoroughly with a spatula and leave it closed container.

**Dough setting and kneading:** Check from time to time whether the material is ready for kneading (if does not stick to gloves and vessel wall). The setting time of the dough is given in Table 1. The setting time is based on an ambient temperature of 20-23 ° C. Higher temperatures shorten and lower temperatures extend the setting time.

When the material is not sticky, it should be thoroughly kneaded within 1-2 minutes. While kneading hands (in gloves) must not be wet and smeared with foreign substances such as cosmetic creams, talc, etc.

**Application onto the model:** Put the acrylic dough on the model and form the orthodontic appliance. If the dough is too dry, moisten it with the previously prepared working solution.

**Polymerization (Villocryl Ortho Mix)**

Polymerization: Place the stone model with formed orthodontic appliance in a pressure polymerizer.

Polymerize in water maintaining the parameters given in Table 1.

**Residual monomer level reduction**

Respect the following for the material: recommended mixing ratio, mixing temperature and time of mixing components and parameters of the polymerization process.

Reducing the residual monomer level can be achieved by extending the polymerization time. Before handing over the orthodontic appliance to the dentist, should be soaked in water at room temperature for 48 hours, this operation can also be performed by the patient after the dentist has handed over the denture.

**Final processing**

Use standard method of treatment. Polishing with **Everall's Pumeks** and **Pasta polerska**.

The orthodontic appliance plate should be thoroughly polished. Scratches are a place of dirt accumulation.

Do not wear the appliance during processing and polishing - permanent changes in the structure of the material may occur, leading to an increase in residual monomer level of the finished prosthesis.

**Corrections and repairs of a damaged orthodontic appliance:** Corrections and repairs of a damaged orthodontic appliance can be performed with Villocryl Ortho/Villocryl Ortho Mix, after moistening the corrected or repaired element with monomer. Self-curing denture repair materials - **Villocryl S** or **Villocryl SP** also can be used.

**Packaging and transport of the orthodontic appliance and the dentist's office:** The finished orthodontic appliance should be protected against mechanical damage and high temperature.

The orthodontic appliance may change color or lose its original physical properties when exposed to chemicals or if it is misused.

During use, the orthodontic appliance should be periodically inspected in dental office.

Report any damage of the denture to the dentist. Repair may be only performed by a dental technician.

If a piece of damaged orthodontic appliance is swallowed, consult a doctor immediately. The material after application is not visible on X-ray pictures.

**Side effects**

**Dental technicians**

Irritation, redness or symptoms of hypersensitivity, including the occurrence of anaphylactic shock, in case of allergy to any of the product ingredients.

**Patients**

Irritation, inflammation and allergic reactions caused by residual monomer washed out of the orthodontic appliance.

Additionally, the inflammation of the oral cavity caused by fungi and bacteria adhering to the orthodontic appliance.

**Storage conditions**

Product and color concentrates are sensitive to temperature and UV radiation. Protect against light and UV radiation.

Store the product components in their original, tightly closed containers, in a dark and well-ventilated room, at a temperature of 5-25 ° C out of the reach of third parties.

**Data necessary to identify the product and the contents of the package**

Reference number, product name, color, package size, batch number and expiry date - imprint on the product packaging.

Package content - see table below.

**Assortment\*****Color shades**

Villocryl Ortho	Villocryl Ortho Mix	Color concentrates
Transparent	Transparent	yellow, orange, red, raspberry, purple, green, light green, blue

**Packaging**

	Powder + liquid + color concentrates (yellow+ red + blue) kit	Powder	Liquid
Villocryl Ortho	500g powder + 250ml liquid + 3x12ml color concentrate	500g, 2 kg, 4 kg	250 ml, 500 ml, 1 l
Villocryl Ortho Mix	500g powder + 250ml liquid + 3x12ml color concentrate	750g, 2 kg, 4 kg	400 ml, 1 l
Color concentrates	3 x 12 ml		50 ml

\* Information on the availability of the assortment is provided by the distributor.

**Application parameters (1)**

Application parameters	Villocryl Ortho	Villocryl Ortho Mix
1 Mixing ratio powder/liquid	-	24 g / 10 g
2 Setting time (23±2 °C)	-	6 - 7 minutes
3 Suitability time for use	-	ok. 15 minutes
4 Polymerization	Min. 20' - 50-60°C 2 atm	Min. 20' - 45-55°C 2 atm

**Problem prevention (2)**

See website: [www.everall.pl](http://www.everall.pl)

Please send any comments regarding our products to our address.

Information provided in any way, even during the demonstration, does not invalidate the instructions for use. Users are responsible for checking that the product is suitable for the intended use. The manufacturer is not liable for damages, including to third parties, resulting from non-compliance with the instructions or non-use. The manufacturer's liability is in any case limited to the value of the products delivered.

**MD**

The symbol means a medical device as defined in EU Regulation 2017/745.

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