**Villacryl Opaker universal metal masking material**  
**Product description**  
**Product characteristics:** acrylic material for masking metal components of removable partial dentures, crowns and dental bridges veneered with acrylic resin.  
**Indications for use:** removable partial dentures, acrylic-veneered crowns and bridges.  
**Clinical benefits:** Restoration of chewing and biting functions impaired by partial loss of dentition.  
**Classification according to EN ISO 10477:** Type 2, Class 1  
 Target group(s) of patients and intended users of the product  
**Users:** Dental technology professionals: dental technicians making dentures, crowns and bridges.

**Patients:** All patients regardless of age and sex with loss of dentition that allows the use of acrylic removable partial dentures, crowns or bridges.  
**Contraindications:** Do not use in case of hypersensitivity to any of the ingredients of the preparation. Use is not recommended in pregnant or breastfeeding women or by dental technicians who are pregnant or breastfeeding.

**Adverse effects**  
 Dental technicians  
 Irritation, redness or hypersensitivity symptoms including anaphylactic shock if allergic to any of the ingredients of the product.  
**Patients**  
 At the moment, no adverse effects have been identified.  
**Product form:** coloured powder in a plastic container, liquid in a glass bottle.  
**Delivery forms:** kit (powder, liquid, etchant) or sold separately: liquid, etchant.

**Product composition**  
**Powder:** acrylic polymer, benzoyl peroxide, pigments  
**Liquid:** acrylic monomers  
**Primer:** silane, isopropyl alcohol  
**Compatible products**  
 For use with metal alloys used in prosthetics and acrylic materials for dentures, dental crowns and bridges.  
**Safety instructions**  
**Risks**  
 H225 Highly flammable liquid and vapour.  
 H317 May cause an allergic skin reaction.  
 H319 Causes serious eye irritation.  
 H336 May cause drowsiness or dizziness.  
 EUH 208 Contient du peroxyde de benzoyle. May produce an allergic reaction.

**Precautionary measures**  
 P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.  
 P261 Avoid breathing dust/fume/gas/mist/vapours/spray.  
 P271 Use only outdoors or in a well-ventilated area.  
 P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection/...  
 P302+P352 IF ON SKIN: Wash with plenty of water/...  
 P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.  
 P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. In case of accidental ingestion, contact a doctor/physician. Do not discard into drains.

**Storage conditions**  
 Store product components in original, tightly closed packaging; avoid direct sunlight; in a well-ventilated room, at 5–25°C, out of the reach of third parties.  
**Waste disposal**  
 Dispose of in accordance with current business waste management regulations.  
 Cardboard packaging, labels, instructions: recyclable, dispose of as non-hazardous separated municipal waste – paper.  
 Ready-to-use disposable wipes soaked in a disinfectant solution can also be used. In addition, the denture should be washed with warm water and soap. Communicate the relevant information to the dentist.  
 Non-hazardous powder: mixed municipal waste Powder packaging: non-hazardous separated municipal waste – plastic. Liquid and liquid packaging – hazardous waste. Hand over to an authorised waste collector.  
 Primer and primer packaging – hazardous waste. Hand over to an authorised waste collector.

**Directions for use**  
**Preparation:** Clean the metal surface of any residual protective mass and metal oxides with a dental sandblaster. Prepare the cleaned metal surface by sandblasting with pure aluminium oxide with a granulation of 250 µm at a pressure of 4 bar. Clean and degrease the prepared surface with acetone, apply Opaker Primer and spread evenly with compressed air (do not use oil compressors). Opaker Primer is used for the acrylic-to-metal chemical bond and prevents the formation of the so-called marginal gap.  
**Mixing:** Pour a few drops of the Opaker Liquid monomer into a separate container and add Opaker Powder until saturation. Mix for approximately 30 seconds until a homogeneous consistency is achieved. Opaker Powders can be mixed together to achieve the desired colour.

**Application:** Coat the metal surface with a thin, opaque layer of the prepared agent using a brush and wait approximately 5 minutes.

**Curing:** The curing of the opaquing agent can take place in a hot air stream or over a burner flame. Care must be taken not to burn the agent's layer. Otherwise, the material will change its structure and will fail to adhere to the metal substrate.  
 Prolonged exposure to temperatures above 250°C may cause the polymer to decompose. Flammable and toxic substances may be released during decomposition.

When properly cured, the agent is hard, has a uniform colour, shows proper adhesion to the metal and masks it accurately.  
**Finishing and polishing:** After curing, wash the surface of Opaker with alcohol to remove the unpolymerised layer.  
**Packaging and transporting dentures/crowns to the dental surgery:** The finished denture, crown or bridge made with the product should be protected from mechanical damage and high temperatures and handed over to the dentist.

**Hygiene, sterilisation, disinfection**  
**Product**  
 The product is designed for reuse. It is not supplied sterile, nor is it required to be sterilised or disinfected before use. Apply the principles of good laboratory practice in dental laboratories to protect the product from contamination.  
**Denture, framework / crown metal component masked with the product**  
 No sterilisation or disinfection procedures are required before the actual acrylic layer of the restoration is placed.  
**Crown / prosthesis made using the product**  
 There is a risk of cross-infection. The crown should be washed and disinfected by the dentist before being placed in the patient's mouth for the first time.

For disinfecting crowns / dentures made of the product, it is recommended to use disinfectants intended by their manufacturers for disinfecting dental prostheses and to follow the manufacturer's recommended procedures.

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**Packaging contents** – see table below.

REF	Set
V210Z10	Villacryl Opaker 4x7 g powder (grey, yellow, brown, pink) + 2x12 ml liquid + 12 ml primer
V210Z12	Villacryl Opaker 2x7 g powder (pink) + 12 ml liquid + 12 ml primer
V210L03	Villacryl Opaker Liquid 12 ml
V210L04	Villacryl Opaker Primer 12 ml

**Prevention of problems (2)**  
 See website [www.everall7.pl](http://www.everall7.pl).

For details on all Everall7 materials, please call: +48 602 701 001.  
 Please send any comments on our products to our address.  
 Information provided in any way, including demonstrations, does not replace reading the instructions for use. Users are required to make sure that the product is suitable for the intended use. The manufacturer shall not be held liable for damage, including to third parties, resulting from a failure to follow the instructions or unsuitability for use. The manufacturer's liability is in any case limited to the value of the products supplied.

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