VILLACRYL SP	INSTRUCTIONS FOR USE - EN
---------------------	----------------------------------

I. Product description

Villacryl SP Acrylic resin for making acrylic parts in framework dentures and complete acrylic dentures with pouring technique.

ISO classification 20795-1: Type 2, Class 2.

Processing parameters: necessary parameters for the proper processing of the materials are given in **Table 1**.

Before reusing the product, check if the fluid retains its properties. In the event of visible traces of density change, cloudiness or discoloration, stop using the product.

Preparation:

Teeth preparation

Use classic acrylic teeth. During setting, the teeth should have a grinded mucosal surface and retention holes / grooves made with a dental milling cutter. Thoroughly clean the mucosal surface of the acrylic teeth from the **WAX RESIDUE** and moisten it with monomer. In special cases, follow the recommendations of the tooth manufacturer.

Product preparation

Shake the powder and liquid containers 2 or 3 times as shown on the package and product labels to help disperse the powder in the monomer. Open the powder and liquid containers 2 times on the packaging and product labels.

I. Making acrylic parts of frameworks.

In order to make the metal parts of the partial denture, it is recommended to use **Villacryl Opaker**.

Preparation: Set the teeth in the wax and make the mask using technical putty silicone. Remove mask from the model and scald the wax. Fix the teeth in the mask. Use a brush to cover the plaster model twice with **Isa-sol** – isolating liquid. Wait for the isolation membrane to form. Cut wide slots (pouring channels) from the palatal or vestibular side of the mask.

Product form

Separate: a mixture of powder in plastic containers and liquid in bottles. The product is available in the form of a set: powder + liquid or separately as a powder and a liquid.

Chemical composition

Powder: polymethyl methacrylate (PMMA), barbituric acid derivative, pigments, viscose fibers *

Liquid: methyl methacrylate (MMA), cross-linking agent, amine catalyst, stabilizer

* Depending on the color variety, the composition of the product may slightly differ in terms of quantity and quality

Product compatibility

For use with acrylic artificial teeth and metal alloys used in dental technology industry

Safety instructions

Warnings

H235 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

The ingredients of the product and the dust generated during the fabrication of dentures may cause mechanical eye irritation. The liquid is sensitive to temperature and UV radiation. An uncontrolled polymerization reaction may occur.

Residual monomer level reduction:

Respect the following for the material: recommended mixing ratio, mixing temperature and time of mixing components and parameters of the polymerization process.

Reducing the residual monomer level can be achieved by extending the polymerization time. Before handing over the denture to the dentist, prosthesis should be soaked in water at room temperature for 48 hours, this operation can also be performed by the patient after the dentist has handed over the denture.

Final processing: Remove the mask and the denture from the model. Use standard method of treatment.

Polishing with Everall’s Pumeks and Pasta polska.

Corrections and repairs of a damaged denture: Corrections and repairs of a damaged denture can be performed with Villacryl SP or self-curing acrylic material - Villacryl 5 after moistening the corrected or repaired elemet with monomer.

II. Making acrylic dentures with pouring technique in agar

Preparation: After modelling the denture on the stone model, place it in water in order to saturate the model (with water), then place it in a special agar flask and fix it with wax to the flask’s base. Leave the pouring inlets closed. Close the flask and then fill it with prepared agar by pouring this stream of earlier prepared agar. Cool the flask to room temperature. Open the flask’s base and gently remove the model to avoid damaging the agar mold. Remove the teeth from the wax denture and scald it to degrease and dean the teeth from wax residues. Scald the stone model and cover the plaster model **twice with Isa-sol** – isolating liquid with a brush. Wait for isolation membrane to form. In the agar mold, cut out the pouring channels for the acrylic. Leave the pouring inlets open. After the model has cooled down, rinse the teeth with the monomer and place them in their places in the agar mold. Place the model on the flask’s base and close it. Place the flask with the pouring inlets facing up.

Measuring and mixing

Recommended mixing ratio: (Table 2) It is essential to follow the stated mixing ratio of the components.

Not following the instructions may result in the formation of a porous structure, a reduction in mechanical properties and increase the level of residual monomer in the finished prosthesis.

Mixing: Pour the appropriate amount of liquid into a chemically resistant vessel and carefully measure the amount of powder. Mix powder and liquid thoroughly with a spatula.

Pouring: After mixing, wait until the consistency of liquid honey is obtained, then pour into the flask a thin stream of acrylic resin. Wait until a matt skin forms on the acrylic surface and polymerize.

Polymerization: Use pressure polymerization in polymerization pot filled with water.

The product is not intended for use with techniques such as, for example, curing with light, microwave or injection method.

The correct parameters of the polymerization process for the material are given in **Table 1**.

Cooling and deflasking

The correct parameters of the cooling process after polymerization for the material are given in **Table 1**.

They guarantee minimal shrinkage after polymerization, which has an impact on the dimension stability of the prosthesis and its retention to the palate.

Final processing: Remove the acrylic pouring canal. Use standard method of treatment. Polishing with **Everall’s Pumeks and Pasta polska**.

The denture plate should be thoroughly polished. Scratches are a place of dirt accumulation. Do not overheat the material during processing and polishing – permanent changes in the structure of the material may occur, leading to an increase in residual monomer level in the finished prosthesis.

Residual monomer level reduction:

Respect the following for the material: recommended mixing ratio, mixing temperature and time of mixing components and parameters of the polymerization process.

Reducing the residual monomer level can be achieved by extending the polymerization time. Before handing over the denture to the dentist, prosthesis should be soaked in water at room temperature for 48 hours.

Corrections and repairs of a damaged denture: Corrections and repairs of a damaged denture can be performed with Villacryl SP or self-curing acrylic material - Villacryl 5 after moistening the corrected or repaired element with monomer.

Packaging and transport of the denture to the dentist’s office: The finished denture should be protected against mechanical damage and high temperature and handed over to the dentist.

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Produkt

The product is intended for multiple use. It is not supplied sterile and does not require sterilization or disinfection prior to use. To protect the product against contamination, apply the principles of good laboratory practices in dental laboratories.

The prosthesis made of materials.

There is a risk of transmission of infection. The prosthesis should be cleaned and disinfected by the dentist before being placed in the patient’s mouth for the first time.

For the disinfection of dentures made of materials, it is recommended to use disinfectants for disinfecting dental prostheses, following the procedures recommended by the manufacturer.

It is also possible to use agents based on at least 65% ethyl or propyl alcohol. The denture should be sprayed with a disinfectant solution, and then gently wiped dry with a paper tissue. You can also use ready-to-use disposable wipes soaked in a disinfecting solution. Additionally, the denture should be washed with warm soapy water. Provide the appropriate information to the dentist.

Batch number and expiry date – see imprint on product labels and packaging.

The expiry date refers to properly stored unopened packaging. Do not use the product after the expiry date.

Additional information

Responsible of setting incidents: Any serious incidents related to the products used do not constitute must be reported to the manufacturer and the competent authority in accordance with local regulations.

Data necessary to identify the product and the contents of the package

Reference number, product name, color, package size, batch number and expiry date – imprint on the product packaging.

Package content – see table below.

Assortment*

Color shades

Villacryl SP

0 - transparent

V2 – milk pink veined

V4 – pink veined

Packaging

	Powder + liquid kit	Powder	Liquid
Villacryl SP	500g powder + 300 ml liquid;	500 g	300 ml

* Information on the availability of the assortment is provided by the distributor.

Application parameters	Frameworks	Pouring technique in agar
1 Mixing ratio powder/liquid	10 g / 6,7 g (7 ml)	10 g / 5 g (5,2 ml)
2 Mixing time (23±2°C)	30” - 60”	30” - 60”
3 Polymerization	min. 20” 50°C - 60°C Pressure 2 bar	30” 60°C Pressure 2 bar
4 Cooling	Remove the model from the polymerizer and cool to room temperature.	Remove the flask from the polymerizer and cool to room temperature.

Problem prevention (2)

See website: **www.everall.pl**

Application parameters	Frameworks	Pouring technique in agar
1 Mixing ratio powder/liquid	10 g / 6,7 g (7 ml)	10 g / 5 g (5,2 ml)
2 Mixing time (23±2°C)	30” - 60”	30” - 60”
3 Polymerization	min. 20” 50°C - 60°C Pressure 2 bar	30” 60°C Pressure 2 bar
4 Cooling	Remove the model from the polymerizer and cool to room temperature.	Remove the flask from the polymerizer and cool to room temperature.

* Information on the availability of the assortment is provided by the distributor.

Application parameters	Frameworks	Pouring technique in agar
1 Mixing ratio powder/liquid	10 g / 6,7 g (7 ml)	10 g / 5 g (5,2 ml)
2 Mixing time (23±2°C)	30” - 60”	30” - 60”
3 Polymerization	min. 20” 50°C - 60°C Pressure 2 bar	30” 60°C Pressure 2 bar
4 Cooling	Remove the model from the polymerizer and cool to room temperature.	Remove the flask from the polymerizer and cool to room temperature.

Problem prevention (2)

See website: **www.everall.pl**

Application parameters	Frameworks	Pouring technique in agar
1 Mixing ratio powder/liquid	10 g / 6,7 g (7 ml)	10 g / 5 g (5,2 ml)
2 Mixing time (23±2°C)	30” - 60”	30” - 60”
3 Polymerization	min. 20” 50°C - 60°C Pressure 2 bar	30” 60°C Pressure 2 bar
4 Cooling	Remove the model from the polymerizer and cool to room temperature.	Remove the flask from the polymerizer and cool to room temperature.

Problem prevention (2)

See website: **www.everall.pl**

Application parameters	Frameworks	Pouring technique in agar
1 Mixing ratio powder/liquid	10 g / 6,7 g (7 ml)	10 g / 5 g (5,2 ml)
2 Mixing time (23±2°C)	30” - 60”	30” - 60”
3 Polymerization	min. 20” 50°C - 60°C Pressure 2 bar	30” 60°C Pressure 2 bar
4 Cooling	Remove the model from the polymerizer and cool to room temperature.	Remove the flask from the polymerizer and cool to room temperature.

Problem prevention (2)

See website: **www.everall.pl**

Application parameters	Frameworks	Pouring technique in agar
1 Mixing ratio powder/liquid	10 g / 6,7 g (7 ml)	10 g / 5 g (5,2 ml)
2 Mixing time (23±2°C)	30” - 60”	30” - 60”
3 Polymerization	min. 20” 50°C - 60°C Pressure 2 bar	30” 60°C Pressure 2 bar
4 Cooling	Remove the model from the polymerizer and cool to room temperature.	Remove the flask from the polymerizer and cool to room temperature.

Information provided in any way, even during the demonstration, does not invalidate the instructions for use. Users are responsible for checking that the product is suitable for the intended use. The manufacturer is not liable for damages, including to third parties, resulting from non-compliance with the instructions or non-use. The manufacturer’s liability is in any case limited to the value of the products delivered.

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