VILLACRYL THERMO PRESS INSTRUCTIONS FOR USE - EN

Thermoplastic material for denture bases

Product characteristics: Thermoplastic material in the form of a granulated product for making dentures using the pressure injection method, intended for use in dental and technical laboratories. The product does not contain methyl methacrylate.

Classification according to EN ISO 20795-1: Typ 3

Chemical composition

Granulated product: mikrocrystaline polyamide, pigments

Performance: Technical data necessary for correct material processing are given in Table 1.

Intended use of the product

Indications: total dentures (edentulous), partial, retention elements in framework dentures, micro dentures, splints and space maintainers.

Contraindications: Do not use the product in patients who are allergic to any of the ingredients.

Target group or groups of target patients and intended users: People of different ages, with partial or complete lack of teeth, eligible for the prosthetic treatment with removable artificial teeth.

Side effects: Not found in the current period of the use of the product. In the absence of proper hygiene, irritations and oral inflammations caused by fungi and bacteria adjacent to the denture are possible.

Clinical benefits: Restoring the chewing function, correct phonetics and improving facial aesthetics, which positively affects the patient's psychological comfort.

Directions for use

Preparation of the working model: The working model should be made of the IV grade dental plaster, e.g. **Stodent IV**, which will ensure perfect reproduction and fitting in the patient's mouth. Before starting the design of the denture, the paralelometer analysis should be performed, in order to determine the insertion path and examine the depth of the undercuts on the model. Thanks to the flexibility of the material, it is possible to perform works using the undercuts of the gum tissue and the abutment tooth for retention. When designing the retention elements in the front section, minimal tooth coverage is recommended for the sake of aesthetics. In the case of lateral teeth, the use of larger clasps will allow for greater stability of the denture. In order to achieve satisfactory retention, the undercuts of the gum tissue and teeth should be used. It is most convenient to use a pencil to apply the design to the model.

The material injection procedure recommended by the manufacturer of the device in which the material will be injected must be strictly observed!
Note! Place the flask in the injection molding machine at a temperature from 90°C to 100°C.
Note! Do not set the injection temperature higher than 285°C.
Unflasking: After the injection process, let the flask cool down for about 30 minutes. Do not speed up the cooling process with water or other coolants.

Note! Do not hit the flask, do not damage its walls during unflasking. The damaged flask is not suitable for injections and cannot be used safely. Any cracks, scratches, deformations on the flask walls or bending the flask axis means that it should be replaced with a new one. Performing work using the defective flask threatens with the risk of health damage.

Pressing: Using the separating disk (e.g. openwork with double-sided diamond embankment), cut off the injection channels. Use the standard method of processing for thermoplastics 45hine TC CUTTERS.

Note! It is not recommended to process the material at high speed. Avoid overheating the material. **Preliminary finishing:** For initial processing, it is recommended to use burs designed for thermoplastic materials. Perform preliminary polishing using rubber polishes and pumice (**Pumice** series by Everall7 and the **45hine Polishing Powder Hard line**), applying low polishing speeds and continuous moistening of the denture to prevent overheating.

Final polishing: For final polishing, it is recommended to use Everall7 Polishing paste / 45hine Polishing Paste Thermoplastic or pastes dedicated to thermoplastic materials. Perform the procedure dry, at low speeds, to prevent overheating of the denture.

Caution should be exercised during polishing: use protective clothing, safety goggles, and dust masks when working with the polishing unit. Use a local extraction system with appropriate power and dust filtration.

Packaging and transport of the denture to the dentist's office: The prepared denture should be protected against mechanical damage and high temperature.

Corrections and repairs: The Villacryl Thermo Press does not chemically bond with traditional acrylic materials. Corrections and repairs of the denture made of the Villacryl Thermo Press may be performed by material injection with the use of fluids for repairs and relinings of polyamide-based thermoplastics. In order to perform the material injection in the case of repair or relining, apply a small amount of material for polyamide-based repairs before performing the injection to the previously developed edges of the denture.

Additional notes: The dentures should be polished thoroughly. The scratches left behind are the place where dirt accumulates. Avoid modeling the denture plate in such a way that the thermoplastic material becomes increasingly thinner, below a thickness of 1 mm.

Reprocessing the material is prohibited. Do not overheat the material during processing and polishing - there may be permanent changes in the material structure.

Do not mix this thermoplastic material with other materials.

In case of any allergic reaction, stop working immediately and consult the physician.

Heats and counteraction

Information for dental technicians: In the case of skin contact with the material in a liquid state, while melting, there is a risk of skin burns. There is the possibility of formation of dangerous vapors after exceeding the recommended processing temperature of the product.

In the case of contact with gases released while melting and dust generated during processing, there is the possibility of irritating the eyes and airways.

In the case of inhalation of gases released in the case of product decomposition, there is the possibility of irreversible effects caused by inhalation.

In the case of inhalation or contact with dust generated while processing, there is the possibility of breathing disorders and skin lesions.

the end. In most devices, the cartridge is inserted with the cap facing upwards (for crimping systems).

In some devices with a plunger system, a brass piston is additionally placed on the cap side, and the injection is performed from the longer end of the tube.

NOTE! It is not recommended to fill the tube with too large an amount of material (max. 25 g) due to the limited heating range/height of the heating elements in injection units.

Material injection: Exercise particular caution during material injection. Use protective clothing, gloves, and safety goggles when operating the injection unit. Failure to follow these recommendations may result in burns.

Injection parameters for Villacryl Thermo Press: The following parameters should be observed to achieve optimal properties of the finished denture. However, they may vary depending on the performance and type of injection unit used.

NOTE! Remember to place the prepared cartridge in the device only after the set temperature has been reached. Do not heat the cartridge from a cold furnace, as this prolongs the total heating time and may cause changes in colour and brittleness of the injected denture (applies mainly to the J100 Pressing Dental unit – for the Polish version only).

Table 1.

Technical parameters for injection	
Injection temperature [°C]	275 - 285
Melting time [min]	15 - 27
Small cartridge	275 °C 15 min
Medium cartridge	275 °C 19 min
Big cartridge	285 °C 24 min
Injection pressure [bar]	5 – 9.0
Cooling time [min]	min. 15

The material injection procedure recommended by the manufacturer of the device in which the material will be injected must be strictly observed!

Note! Place the flask in the injection molding machine at a temperature from 90°C to 100°C.

Note! Do not set the injection temperature higher than 285°C.

Unflasking: After the injection process, let the flask cool down for about 30 minutes. Do not speed up the cooling process with water or other coolants.

Note! Do not hit the flask, do not damage its walls during unflasking. The damaged flask is not suitable for injections and cannot be used safely. Any cracks, scratches, deformations on the flask walls or bending the flask axis means that it should be replaced with a new one. Performing work using the defective flask threatens with the risk of health damage.

Pressing: Using the separating disk (e.g. openwork with double-sided diamond embankment), cut off the injection channels. Use the standard method of processing for thermoplastics 45hine TC CUTTERS.

Note! It is not recommended to process the material at high speed. Avoid overheating the material. **Preliminary finishing:** For initial processing, it is recommended to use burs designed for thermoplastic materials. Perform preliminary polishing using rubber polishes and pumice (**Pumice** series by Everall7 and the **45hine Polishing Powder Hard line**), applying low polishing speeds and continuous moistening of the denture to prevent overheating.

Final polishing: For final polishing, it is recommended to use Everall7 Polishing paste / 45hine Polishing Paste Thermoplastic or pastes dedicated to thermoplastic materials. Perform the procedure dry, at low speeds, to prevent overheating of the denture.

Caution should be exercised during polishing: use protective clothing, safety goggles, and dust masks when working with the polishing unit. Use a local extraction system with appropriate power and dust filtration.

Packaging and transport of the denture to the dentist's office: The prepared denture should be protected against mechanical damage and high temperature.

Corrections and repairs: The Villacryl Thermo Press does not chemically bond with traditional acrylic materials. Corrections and repairs of the denture made of the Villacryl Thermo Press may be performed by material injection with the use of fluids for repairs and relinings of polyamide-based thermoplastics. In order to perform the material injection in the case of repair or relining, apply a small amount of material for polyamide-based repairs before performing the injection to the previously developed edges of the denture.

Additional notes: The dentures should be polished thoroughly. The scratches left behind are the place where dirt accumulates. Avoid modeling the denture plate in such a way that the thermoplastic material becomes increasingly thinner, below a thickness of 1 mm.

Reprocessing the material is prohibited. Do not overheat the material during processing and polishing - there may be permanent changes in the material structure.

Do not mix this thermoplastic material with other materials.

In case of any allergic reaction, stop working immediately and consult the physician.

Heats and counteraction

Information for dental technicians: In the case of skin contact with the material in a liquid state, while melting, there is a risk of skin burns. There is the possibility of formation of dangerous vapors after exceeding the recommended processing temperature of the product.

In the case of contact with gases released while melting and dust generated during processing, there is the possibility of irritating the eyes and airways.

In the case of inhalation of gases released in the case of product decomposition, there is the possibility of irreversible effects caused by inhalation.

In the case of inhalation or contact with dust generated while processing, there is the possibility of breathing disorders and skin lesions.

Information for patients and dental technicians:

The material is not resistant to chemicals with pH<4.

The material is not resistant to long-term effect of alcoholic beverages with concentrations above 40%.

The material is not resistant to oxidizing chemicals, e.g. preparations, such as COREGA TABS.

The material is not visible on X-ray.

Information for patients:

Information for patients to be provided by the dentist after the denture is made:

The denture is not resistant to oxidizing chemicals, e.g. preparations, such as COREGA TABS, thus it is not allowed to use cleaning agents, i.e. COREGA TABS (or its equivalents), as well as mouthwashes.

Wash the denture after each meal with a soft brush under running water. For cleaning, you can use pastes for standard hygiene of dentures that do not contain abrasive substances.

It is recommended to wash the dentures from time to time in the ultrasonic cleaner here filled with a special cleaning liquid for dentures made of thermoplastics.

Wash the unused denture and store it in a sealed container. It is not allowed to leave the denture in the sun and near heat sources. It is forbidden to place the denture in an open fire - flammable material.

The denture is not resistant to long-term effect of alcoholic beverages with concentrations above 40%. The denture is not visible on X-ray.

The denture can change color or lose its original physical properties in the case of contact with chemicals or improper use.

In the case of any allergic reaction to the material from which the denture is made, immediately discontinue the use and consult the physician.

When using it, the patient should provide the denture for a periodic dental check.

It is recommended to reline the denture if it lacks proper retention due to the loss of the alveolar process.

Information for dentists:

It is possible to transmit the infection. Each denture should be washed and disinfected before placing in the mouth. It is advisable to inform the user of the denture about the need to apply hygiene rules related to the use of the denture. For dentures' hygiene, special brushes for cleaning dentures, unscented soap, preparations intended for cleaning and care of dentures (except for preparations with oxidizing properties) and ultrasonic washers should be used. Do not use abrasives.

Storage and handling

The product is sensitive to UV radiation. Store the product in the original airtight packaging and in dark and well-ventilated room at a temperature from 5°C to 25°C, as well as in the place inaccessible to third parties. Open the product immediately before injection. Do not use the product after the expiration date stated on the packaging and labels of the product ingredients. Before using the product, check whether the bag with moisture absorber is damaged.

Waste management

Cardboard packaging, labels, instructions - the possibility of recycling, dispose of as non-hazardous municipal waste.

The granulated product and the granulated product packaging are not hazardous waste. Dispose of as non-hazardous municipal waste. In the case of the packaging, the possibility of recycling.

Residues of the material after injection should be treated as municipal waste.

Note!

Material for professional use only in technical and dental laboratories.

Information provided in any way, even during demonstrations, does not violate the instructions for use. The users are obliged to check whether the product is suitable for the intended use. The manufacturer is not liable for damages, including third parties, resulting from non-observance of the instructions or unsuitability for the application. The manufacturer's responsibility is in every case limited to the value of the delivered products.

Please send any comments concerning our products to our address.

Each serious incident concerning the product should be reported to the manufacturer and the competent authority of the country where the user or the patient resides.

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