VILLACRYL DISC DB/CLEAR INSTRUCTIONS FOR USE - EN

Villacryl DISC DB – PMMA disc for milling removable denture bases

Villacryl DISC CLEAR - PMMA discs for orthodontic applications

Product description

Villacryl DISC DB are PMMA discs designed for CAD/CAM technology to produce removable denture plates.

Villacryl DISC CLEAR are PMMA discs designed for CAD/CAM techniques, used to make occlusal splints, relaxation splints and surgical templates.

Indications for use: Complete dentures (edentulism), partial removable dentures.

Occlusal splints, relaxation splints and surgical templates.

Clinical benefits: Restoration of chewing and biting functions lost due to partial or total loss of natural teeth; orthodontic treatment.

Target group(s) of patients and intended users of the products:

Users: Professionals in the dental technology industry - certified dental technicians making dental prostheses using CAD-CAM technology.

Patients: people of all ages and genders, with partial or complete loss of teeth, eligible for prosthetic treatment whereby removable dentures and orthodontic treatment including children, adolescents, and adults.

Contraindications

Do not use the product for people with a diagnosed allergy to the ingredients of the product. It has side effects, use a different material, suitable for the patient.

Adverse effects

Dental technicians

Irritation, redness or hypersensitivity symptoms including anaphylactic shock if allergic to any of the ingredients of the product.

Patients

Irritation, inflammation and allergic reactions of the oral cavity caused by fungi and bacteria adhering to the denture.

Product form

Villacryl DISC DB: Disc 25 or 30 mm thick and 98.5 mm in diameter, either clear or in various hues of pink.

Villacryl DISC CLEAR: 16 and 20 mm thick transparent disc.

Forms of delivery: PMMA disc wrapped in a cardboard box and protective foam.

Product composition: methyl methacrylate polymer, pigments. Some colour variations may contain viscose fibres.

Compatible products

For use with acrylic artificial teeth and materials for bonding artificial teeth into a milled denture plate.

Safety instructions

Risks

When coming into contact with the technician's body, dust generated during milling can irritate the skin, cause an allergic skin reaction or irritate the respiratory tract.

If dust particles enter the eyes, mechanical irritation of the eyes may result

Precautions

Avoid inhaling dust generated during milling.

Wear protective gloves/protective clothing/eye protection/face protection.

Do not eat or smoke while working with the material.

Use ventilation in the workplace.

If dust comes into contact with eyes, flush using water and consult a doctor.

Wash areas of skin that have come into contact with the product ingredient thoroughly using soap and water.

Information to be provided to patients / denture users:

If symptoms of allergic reactions or irritation occur, discontinue use of the denture or orthodontic splints immediately and consult a doctor.

Strong oxidizing agents contained in cleaning products and mouthwashes may change the color of the dentation during use.

After each meal, the denture or an orthodontic splint should be washed with a soft toothbrush under

cool running water. Unscented soap, various denture cleaners and denture care products that do not contain abrasive substances can be used for cleaning.

We recommend to wash the dentures or an orthodontic splint periodically in an ultrasonic cleaner filled with denture cleaning liquid. When not used, wash and store denture in a special container designed for the storage of dentures.

Do not leave the denture or an orthodontic splint under direct sunlight or near heat or open flame.

The prosthesis is not resistant to prolonged exposure to alcoholic beverages.

A prosthetic or an orthodontic splint restoration may change colour or lose its original physical properties if it comes into contact with chemicals or is misused.

During use, the patient should attend periodic check-ups of the denture by the dentist.

Any damage to the denture should be reported to the dentist. Repairs can only be carried out by a dental technician.

If a piece of damaged denture or an orthodontic splint is swallowed, see a doctor immediately. The material from which the restoration is made is not visible in X-rays. In case of atrophy of the alveolar process, we recommend to pad or rebase the dental prosthesis plate.

Storage conditions

Store the product in the original packaging in a dark room, at a temperature of 5-25 °C, away from unauthorized persons.

Waste disposal

Dispose of in accordance with current business waste management regulations.

Cardboard packaging, labels, instructions: recyclable, dispose of as non-hazardous separated municipal waste - paper.

Acrylic disc: non-hazardous segregated municipal waste - plastic.

Directions for use (step-by-step instructions)

For proper product processing, dental design skills using CAD/CAM are required. The discs should be milled using a compatible CAD/CAM system.

It is important not to overheat the material!

Preliminary activities:

Read the instructions carefully before using the product.

Before making the prosthesis, ensure that the patient for whom the prosthesis or an orthodontic splint is being made is not allergic to any of the material's ingredients.

NOTE: In order to choose the right colour for the denture plate, it is recommended to use the colour matching chart (colour chart). Colour should be chosen individually and if in doubt, do your own test.

Before using the disc, read the disc label carefully to select the correct disc height and hue. To design the prosthesis plate, use appropriate CAD software. When designing, the minimum wall (prosthesis plate) thickness must be 2 mm. Failure to observe the minimum thickness may result in damage (e.g. cracking). Villacryl DISC DB Villacryl DISC CLEAR can be milled in CAD/CAM milling systems using standard 98.5 mm chuck. Only the milling tools recommended for this system should be used to process the material. Note that the material may overheat during dry milling. This can damage the material. In such case, we recommend wet milling.

Bonding teeth to the denture plate

Artificial teeth should be bonded to the denture plate using PMMA self-curing resin (Villacryl S) or photocurable materials, according to the manufacturer's instructions. The mucosal surface in artificial teeth should be roughened and prepared according to the manufacturer's instructions. In addition, bond strength can be increased by applying Villacryl Bond material to the tooth surface and milled denture plate.

Finishing and polishing of the prosthesis/split

After milling, separate the restoration from the disc connectors using a suitable milling cutter. To finish or adjust the shape, use a cross-cut cutter. Particular attention should be paid to the edges of the milled tooth holes and interdental spaces. The required minimum thickness of the prosthesis of 2 mm must be maintained.

Polishing

Polishing is performed as standard, using polishing tools suitable for acrylic polishing materials.

Adjustments and repairs of damaged dentures/ split

Repairs made from Villacryl DISC DB Villacryl DISC CLEAR can be repaired, adjusted or padded using conventional self-polymerising resin for denture plates (e.g. Villacryl S or Villacryl SP), after you first wet the item to be adjusted or repaired using liquid.

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The finished prosthesis/ split should be protected against mechanical damage and high temperatures.

Higiene, sterilisation, disinfection

Product

The product is designed for the fabrication of a single prosthesis. It is not supplied sterile, nor is it required to be sterilised or disinfected before use. Apply the principles of good laboratory practice in dental laboratories to protect the product from contamination.

Prosthesis/split made using the product

There is a risk of cross-infection. The prosthesis/split should be washed and disinfected by the dentist before being placed in the patient's mouth for the first time.

For disinfecting dentures (split) made of the product, it is recommended to use disinfectants intended by their manufacturers for disinfecting dental prostheses and to follow the manufacturer's recommended procedures.

Agents based on at least 65% ethyl or propyl alcohol can be used. The denture (split) should be sprayed with a disinfectant solution and then gently wiped dry with a paper tissue. Ready-to-use disposable wipes soaked in a disinfectant solution can also be used. In addition, the denture should be washed with warm water and soap. Communicate the relevant information to the dentist.

Additional information

Reporting serious incidents: Any serious incidents involving the products must be reported to the manufacturer and to the competent authority in accordance with local regulations.

Data necessary for the identification of the product and package contents

Reference number, product name, colour, disc thickness, batch number and expiry date - printed on the product packaging.

The expiry date refers to the unopened packaging properly stored. Do not use the product after the expiry date.

Packaging contents

VILLACRYL DISC DB

Colour	Available thickness in mm
V4	25, 30
V26	25, 30
V27	25, 30

Colour	Available thickness in mm
transparent	16, 20

Prevention of problems

See website www.everall7.pl.

Please send any comments on our products to our address.

For details on all Everall7 materials, please call: +48 602 701 001

Information provided in any way, including demonstrations, will not invalidate the instructions for use. Users are required to make sure that the product is suitable for the intended use. The manufacturer shall not be held liable for damage, including to third parties, resulting from a failure to follow the instructions or unsuitability for use. The manufacturer's liability is in any case limited to the value of the product supplied.

MD

The symbol means a medical device as defined in EU Regulation 2017/745.

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