1. Product description
Villacryl S: Self-curing acrylic resin for making dentures repairs.
ISO classification 20795-1: Type 2, Class 1.
Processing parameters: necessary parameters for the proper processing of the materials are given in **Table 1**.
Indications for use
Villacryl S: repair of removable complete and partial dentures, acrylic parts of framework dentures and acrylic orthodontic appliances.
Clinical benefits:
Restoring biting and chewing functions of dentures damaged during use.
Restoring the proper functioning of orthodontic appliances damaged during use in order to correct the dentition or prevent movement of corrected teeth.
Target group or target groups of patients and intended uses of the devices
Users: Professionals in the dental technology industry - certified dental technicians producing dental prostheses.
Patients: People of different ages and sexes, denture users. Uses of orthodontic appliances: children, adolescents, adults.
Contraindications
Do not use the products by people diagnosed with a lack of resistance to allergenic factors, including the ingredients of the product. In such situations, use a different material suitable for the patient.
Do not use the product by pregnant women and nursing mothers in the group of dental technicians.

Product form
Separate a mixture of powder in plastic containers and liquid in bottles.
The product is available in the form of a set: powder + liquid or separately as a powder and a liquid.
Chemical composition
Powder: polymethyl methacrylate (PMMA), barbituric acid derivative, pigments, viscose fibers *
Liquid: methyl methacrylate (MMA), cross-linking agent, amine catalyst, hydroquinone
* Depending on the color variety, the composition of the product may slightly differ in terms of quantity and quality
Product compatibility
For use with acrylic artificial teeth and acrylic dentures.
Safety instructions
Warnings
H225 Highly flammable liquid and vapour.
H315 Causes skin irritation.
H317 May cause an allergic skin reaction.
H335 May cause respiratory irritation.
The ingredients of the product and the dust generated during the fabrication of dentures may cause mechanical eye irritation. The liquid is sensitive to temperature and UV radiation. An uncontrolled polymerization reaction may occur.
Precautions
In case of any allergic reaction, stop work immediately and consult a dermatologist.
P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P302+P352 IF ON SKIN: Wash with plenty of water.
P501 Dispose of the contents/container to an authorized consignee of the waste.
It is recommended to use gloves made of HDPE polyethylene for hand protection, a filtering fluid-compatible face mask for respiratory protection, goggles and workwear.
It is advisable to use forced ventilation of rooms and workplace exhaust systems.
In the event of contact of the product ingredient or dust with the eyes, rinse them with water. Contact a doctor is recommended in case of irritation.
Do not eat or drink while working with the material. In the event of accidental ingestion of an ingredient of the product, consult a doctor immediately.
Do not pour the liquid down the sanitation system.
Before reusing the product, check if the fluid retains its properties. In the event of visible traces of density change, cloudiness or discoloration, stop using the product.

Information to be provided to patients / prosthesis users:
To reduce the residual monomer, it is advisable during the first three days, to soak the prosthesis in water when not in use.
If any symptoms of an allergic reaction or irritation occur, discontinue using the prosthesis immediately and consult a doctor.
Strong oxidizing compounds contained in cleaning agents and mouthwashes may change the color of the prosthesis after polymerization. After each meal, the prosthesis should be washed with a soft brush under cool, running water.
For cleaning, can be used unscented soap, various preparations for cleaning and caring for dental prostheses that do not contain abrasives. It is recommended to wash the dentures from time to time in a vibration washbasin filled with a special cleaning liquid for dentures. The material is not resistant to long-term effects of alcohol beverages.
The repair material may change color or lose its original physical properties when exposed to chemicals or if it is misused.
Report any damage of the denture to the dentist. Repairs may be only performed by a dental technician.
If a piece of damaged dentures is swallowed, consult a doctor immediately. The material after polymerization is not visible on X-ray pictures. It is recommended to rely on rebase the denture base in the event of the disappearance of the alveolar process.

Side effects
Dental technicians
Irritation, redness or symptoms of hypersensitivity, including the occurrence of anaphylactic shock, in case of allergy to any of the product ingredients.
Patients
Irritation, inflammation and allergic reactions caused by residual monomer washed out of the repaired denture.
Storage conditions
Product is sensitive to temperature and UV radiation. Protect against light and UV radiation.
Store the product components in their original, tightly closed containers, in a dark and well-ventilated room, at a temperature of 5-25 °C out of the reach of third parties.
Waste disposal
Dispose in accordance with applicable regulations of economic activity waste management.
Cardboard packages, labels, instructions: non-hazardous segregated municipal waste - paper.
Powder: non-hazardous mixed municipal waste.
Powder packaging: non-hazardous segregated municipal waste - plastic.
The liquid and its packaging - hazardous waste, comply with the regulations on hazardous waste management (storage, transport and disposal) and hand it over to a professional hazardous waste consignee.

How to use (Step by step instructions)
Initial steps:
Read the instructions carefully before using the product.
Make sure the patient for whom the prosthesis is made has not been diagnosed with the lack of resistance to allergenic factors. Before reusing the product, check if the fluid retains its properties. In the event of visible traces of density change, cloudiness or discoloration, stop using the product.

Preparation:

Cleaning and disinfecting dentures before repair

Remove any plaque and scale from the surface of the repaired denture. Use protective gloves and mask.

Denture Cleaner by Everall7 is recommended for cleaning the denture. Wash the cleaned denture with warm water.

Preparation of the prosthesis for repair

When adding teeth to a partial denture, use classic acrylic or composite teeth. The teeth during fixing should have a grounded mucosal surface. Thoroughly clean the mucosal surface of the acrylic teeth from the isolator and moisten it with monomer.

Prepare the repaired prosthesis area thoroughly using the standard method.

Clean the denture's damaged area and moisten it with monomer.

Measuring and mixing

Recommended mixing ratio: (Table 1) It is essential to follow the stated mixing ratio of the components. Not following the instructions may result in the formation of a porous structure, a reduction in mechanical properties and increase the level of residual monomer in the finished prosthesis.

Mixing: Pour the appropriate amount of liquid into a chemically resistant vessel and carefully measure the amount of powder. Mix powder and liquid thoroughly with a spatula and leave in closed vessel for 1-2 minutes in 20-23°C.

Application: The setting time of the dough is given in Table 1 and relates to an ambient temperature of 20-23 °C. Higher temperatures shorten and lower temperatures extend the setting time. Check from time to time whether the material is ready, especially when the ambient temperature is above 20-23 °C. If the dough dries too quickly (high temperature), moisten the surface of the dough placed on the denture with the liquid.

Polymerization:

The product is a material polymerizing at ambient temperature.

For safety reasons (lower residual monomer content, better mechanical parameters) it is recommended to perform pressure polymerization in a repair flask or a pressure polymerization under the conditions given in Table 1.

The product is not intended for use with techniques such as, for example, curing with light, microwave or injection methods.

Ambient temperature polymerization: leave the repaired denture at ambient temperature and wait for the complete curing of the material in repair area.

Water polymerization: It is required to use pressure equipment with a minimum working pressure of 2 bar.

Polymerize the material in water with the parameters given in Table 1.

Residual monomer level reduction:

Respect the following for the material: recommended mixing ratio, mixing temperature and time of mixing components and parameters of the polymerization process.

Reducing the residual monomer level can be achieved by extending the polymerization time.

Final processing: Use standard method of treatment. Polishing with Everall's **Pumex** and **Pasta polska**.

The repaired area should be thoroughly polished. Scratches are a place of dirt accumulation. Do not overheat the material during processing and polishing - permanent changes in the structure of the material may occur, leading to an increase in residual monomer level in the finished prosthesis.

Packaging and transport of the denture to the dentist's office: The repaired denture should be protected against mechanical damage and high temperature and handed over to the dentist.

Hygiene, sterilization, disinfection

Precaution:

The product is intended for multiple use. It is not supplied sterile and does not require sterilization or disinfection prior to use. To protect the product against contamination, apply the principles of good laboratory practices in dental laboratories.

The prosthesis made of materials:

There is a risk of transmission of infection. The repaired prosthesis should be cleaned and disinfected by the dentist before being replaced in the patient's mouth.

For the disinfection of dentures made of materials, it is recommended to use disinfectants for disinfecting dental prostheses, following the procedures recommended by the manufacturer.

It is also possible to use agents based on at least 65% ethyl or isopropyl alcohol. The denture should be sprayed with a disinfectant solution, and then gently wiped dry with a paper tissue. You can also use ready-to-use disposable wipes soaked in a disinfecting solution. Additionally, the denture should be washed with warm soapy water. Provide the appropriate information to the dentist.

Batch number and expiry date: see imprint on product labels and packaging.

The expiry date refers to properly stored unopened packaging. Do not use the product after the expiry date.

Additional information

Reporting of serious incidents: Any serious incidents related to the products must be reported to the manufacturer and the competent authority in accordance with local regulations.

Data necessary to identify the product and the contents of the package

Reference number, product name, color, package size, batch number and expiry date - imprint on the product packaging.

Package content - see table below.

Assortment*

Color shades

Packaging

Villacryl S

Powder + liquid kit

Powder

Liquid

100 g powder + 50 ml liquid

100 g powder + 2 x 50 ml liquid

500 g powder + 250 ml liquid

1 kg

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

200 ml, 500 ml

The symbol means a medical device as defined in EU Regulation 2017/745.

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