

EA everall7



## VILLACRYL ACRYLIC STAINS

Acrylic material for prosthesis characterization



Acrylic resins used for the characterization of dentures are specialized materials applied in dental prosthetics to give dentures a more natural appearance. Thanks to them, it is possible to replicate individual features of the oral tissues, such as gum color or discolorations.

They are important because they improve the aesthetics of prosthetic restorations, which in turn affects the patient's psychological comfort and self-confidence. They are mainly used during the characterization of acrylic dentures as well as in their repair and adjustment. With the use of coloring agents, dentures become less noticeable and better matched to the natural appearance of the oral cavity.



# VILLACRYL ACRYLIC STAINS



Villacryl Acrylic Stains are an acrylic material available as a range of multicoloured powders, intended for the individual characterization of dental prostheses. The material is dedicated to dental technicians who fabricate acrylic dentures for patients with specific aesthetic requirements. Thanks to Villacryl Acrylic Stains, it is possible to give an individualized colouration to the mucosal surface of the denture base as well as to the surfaces of acrylic teeth, including in partial (skeletal) dentures and complete dentures.

The product can also be used in other prosthetic restorations based on acrylic materials, wherever precise and individualized characterization is required. The use of low-translucency colours allows for effective masking of internal metal structures within the denture, improving the aesthetics of the final restoration and achieving a more natural appearance.

## CHARACTERISTICS

- Compatibility of Villacryl Acrylic Stains acrylic resin with the acrylic material used for the denture base
- A wide range of shades enables the reproduction of discolorations and colour variations of the gingiva, from dark to light tones, creating an optical depth effect
- Masking of structural elements – more effective concealment of metal components of the denture or show-through areas

## ADVANTAGES

- Good bonding with acrylic
- A wide range of Villacryl Acrylic Stains shades allows for achieving a natural denture appearance
- The possibility of freely mixing colours offers unlimited customization options
- Powders can be used in both heat-curing techniques (with Villacryl H, Villacryl H Rapid monomers) and cold-curing polymerization techniques (with Villacryl S, Villacryl SP monomers)
- Easy application and processing

### Colours:

#### Basic:

- Light pink
- Dark pink
- Veined pin

#### Stains (tinting colours):

- White
- Transparent veined
- Pink
- Yellow
- Brown
- Blue
- Red
- Dark red

### Available packaging:

VS340Z02

Villacryl Acrylic Stains 3 x 75 g + 8 x 10 g

**Everall7 Sp. z o.o.**  
Augustówka 14  
02-981 Warszawa, Polska  
T +48 22 858 82 72

